



Astell are proud to manufacture highly customisable steam generators that are built-to-order. This enables us to meet any challenges that your facility requirements present. Astell steam generators are designed with storage capacity to meet the high-volume, fast-peak, fast-recovery steam requirements of an autoclave's air removal process. This technical data sheet provides you with a range of measurements and specifications from previous STG072 models. Our Sales and Technical teams will provide you with the precise specification of your steam generator after consultation.

## 1. STG072 General Overview

Steam Outlet	1" Gate Valve
Maximum Steam Flow Rate	108 ( kg/hr )
Maximum Steam Pressure	650 ( kPa ) 6.5 ( bar )
Number of Heaters	3 ( x 24kW )


## 2. The Pressure Vessel

The pressure vessel is where the steam is generated.

Internal Volume	217 ( L )
Operational Water Volume	146 ( L )



A typical STG072 with optional blowdown vessel (left) and water softener (right)



A typical STG072 rear view

## You Will Need

The steam generator has a number of requirements that must be catered for prior to installation. These are: an adequate route between the delivery and installation sites, physical space for the installation, the correct connection points for the appropriate quality of water and electricity, and a suitable drainage connection. Read this section carefully and ensure all the requirements can be met prior to ordering your steam generator.

## 3. Space For The Steam Generator...

	Dimensions
Width	802 ( mm )
Height	2000 ( mm ) ( ± 50 mm )
Depth	1021 ( mm )
Footprint / Area	0.82 ( m <sup>2</sup> )
Clearance Space	500 ( mm )
Dry Weight	500 ( kg )

## ...And For Some Options

	Water Softener	Blowdown Vessel
Width	250 ( mm )	345 ( mm )
Height	500 ( mm )	1050 ( mm )
Depth	440 ( mm )	335 ( mm )

It is recommended that a gap of 500mm is left around the steam generator for maintenance and utilities access.

Remember to leave space for any options you select that are physically external to the steam generator, such as a water softener or blowdown vessel.

A route capable of bearing the size and weight of the steam generator must be available between the delivery and installation site. If such a route is unavailable, discuss the construction with Astell.

## 4. Electricity

	Electrical Requirement
Maximum Voltage	415 ( V )
Wattage	72 ( kW )
Frequency	50/60 ( hz )
Standard Amperage	111 ( A )

	Single	Three
AC Phases	N/A	Available

A three phase steam generator will need to be fitted to an electrical isolator box or commando plug and socket. Plugs and sockets not supplied. Contact Astell or your electrician for more information.

## 5. A Clean Water Supply...

Water Type	Mains Water
Maximum Water Hardness	70 (ppm)
Water pH	pH 7

Is your water quality outside the specification given above? Astell can supply an external water softener kit to reduce the hardness and pH of your mains water to acceptable levels. Please see the table to the right for an overview of the water softener's installation requirements.

Water Enters The Steam generator Via	a washing machine inlet hose with 3/4 inch compression fittings
Minimum Water Flow Rate	5 (L / min)

## ...And Possibly A Water Softener Kit

	Water Softener Kit
Unit Width	250 (mm)
Unit Height	500 (mm)
Unit Depth	440 (mm)
Installation Location	Adjacent to drain AND adjacent to steam generator
Water Input Via	Ø 3/4" male isolator valve
Minimum Water Input Pressure	1.5 (bar)
Maximum Water Input Pressure	6 (bar)
Overflow Connection to Drain	Floor drain or standpipe drain capable of receiving Ø35mm AND Ø26mm hoses simultaneously

## 6. Adequate Drainage for Overflow...

Drainpipe Material	Copper
Drainpipe Diameter	Ø 35 (mm)
Maximum Drain Height	0 (mm) above floor
Drainpipe Connections	Between Steam Generator and Main Building Drain Only

Steam generators can eject precipitate left within them using pressurised steam. This process is recommended to keep steam generators functioning optimally. A blowdown vessel captures the pressurised steam and deposits ejected from the steam generator at velocity, decreasing risk.

A water softener has drainage requirements that are additional to those of the steam generator. If you are fitting a water softener, review the water softener drain requirements in section 4.

If your drainpipe is not made of copper, ask Astell if they can build a steam generator designed for the drains you have.

## ...And Possibly For A Blowdown Vessel

	Blowdown Vessel
Unit Width	345 (mm)
Unit Height	1050 (mm)
Unit Depth	335 (mm)
Installation Location	Adjacent to drain AND adjacent to steam generator
Steam Input Via	Ø 15mm Copper Pipe
Steam Vent Size	Ø 54 (mm)
Drainpipe Diameter	Ø 54 (mm)
Maximum Drain Height	0 (mm) above floor
Drainpipe Connections	Between Steam Generator and Main Building Drain Only

## Services

### 7. Delivery

Astell equipment is sold ex works. For your convenience, Astell will arrange a specialist equipment haulier at cost on your behalf. Alternatively, you may elect to arrange haulage from our headquarters; please instruct our sales team if this is the case.

*For customers in the United Kingdom of Great Britain and Northern Ireland, the following service is also available:*

### 8. Planned Preventative Maintenance

Astell Maintenance ensures that your steam generator is regularly inspected, serviced, and maintained with official parts by a qualified Astell Service Engineer. Speak to your sales representative about building a maintenance contract tailored to your steam generator use.

## Standards

### 9. Applicable Standards

Astell steam generators are manufactured in accordance with the following standards:

- PED 2014/68/EU
- ISO9001:2015
- BS2646

(Manufacturing to other standards is available on request)

All information is accurate on data of issue. In line with continuous advancements in design and manufacturing, Astell reserve the right to modify product specifications. E&OE



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.