

Digital refractometer KERN ORM

NEW



Transport and storage case



Rear view, screw-on battery compartment cover

Digital measurement of refraction index for universal application

Features

- The KERN ORM refractometers are accurate and universal maintenance free digital handheld refractometers
- They are characterized by their easy-using and robustness
- The typical and practical design is suitable for a quick and convenient everyday use
- The large, easy-to-read display with integrated temperature display supports the user to reliably determine the measurement
- The integrated automatic temperature compensation (ATC), avoids the manual conversion of the measurement. This allows a quick and efficient usage of the instrument
- Rapid, user-friendly calibration of the refractometer is possible at any time using standard commercial distilled water
- The refractometers from the KERN ORM range are protected to international IP65 protection class, against dust and water splashes. After use, you can rinse the refractometer under running water
- Mean value measurements possible
- The following accessory-parts are included:
 - Prism cover lid
 - Pipette
 - Storage box
 - 1 x AAA battery
 - Screwdriver

Technical data

- Measurement temperature: 0 °C – 40 °C
- Overall dimensions W×D×H 121×58×25 mm
- Net weight approx. 289 g
- Power supply: 1 x AAA (1,5 V)
- Lifetime of the battery: approx. 10.000 measurements
- ATC (Automatic Temperature Compensation)
- Minimum sample volume: 4 drops
- Automatic energy management (AUTO-OFF after 60 seconds)
- Mean value measurement (15 measurements)

! Also available with calibration certificate, see page 109!

STANDARD




Digital refractometer KERN ORM-B · ORM-R · ORM-SU

Scope of application: Basic measurements for Brix and refractive index

The following models are particularly suitable for basic measurement where the result is required in Brix or refractive index. They are used to determine the sugar content in food or for monitoring processes in the industry (coolant monitoring, water-based mixtures). Alternatively the display can be switched to show Brix or the refractive index.

The main scope of applications is:

- Industry: Monitoring of lubricants in machines and quality control
- Food industry: Beverages, fruits and sweets
- Agriculture: Determination of the degree of ripeness of fruit for quality control in harvesting
- Restaurants and large-scale catering establishment



| Model | Scales | Measuring range | Accuracy | Division |
|-----------------|------------------|--------------------|-------------|-----------|
| KERN | | | | |
| ORM 50BM | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |
| ORM 1RS | Brix | 0 – 90 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,5177 nD | ± 0,0003 nD | 0,0001 nD |

Scope of application: Sugar

The following models are particularly suitable for direct measurement of different types of sugar. These are used to determine the content of the respective type of sugar in water-based liquids. It is possible to switch between the four different scales.

The main scope of applications is:

- Food industry: Beverages, fruits and sweets
- Agriculture: Determination of the degree of ripeness of fruit for quality control in harvesting
- Restaurants and large-scale catering establishment



| Model | Scales | Measuring range | Accuracy | Division |
|----------------|------------------|---------------------|-------------|-----------|
| KERN | | | | |
| ORM 1SU | Fructose | 0 – 69 % | ± 0,2 % | 0,1 % |
| | Glucose | 0 – 60 % | ± 0,2 % | 0,1 % |
| | Brix | 0 – 90 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,577 nD % | ± 0,0003 nD | 0,0001 nD |
| ORM 2SU | Lactose | 0 – 17 % | ± 0,2 % | 0,1 % |
| | Maltose | 0 – 16 % | ± 0,2 % | 0,1 % |
| | Dextran | 0 – 11 % | ± 0,2 % | 0,1 % |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |

Digital refractometer KERN ORM-HO · ORM-NA · ORM-SW

Scope of application: Honey

The following model is particularly suitable for the measurement of the water content in honey according to the International Honey Commission (IHC2002) and “degrees Baumé” to determine the relative density of liquids. Alternatively the display can be switched to show Brix or the refractive index.

The main scope of applications is:

- Beekeeping
- Honey production



| Model | Scales | Measuring range | Accuracy | Division |
|----------------|------------------|--------------------|-------------|-----------|
| KERN | | | | |
| ORM 1HO | Brix | 5 – 38 % | ± 0,2 % | 0,1 % |
| | Baumé | 33 – 48 °Bé | ± 0,2 °Bé | 0,1 °Bé |
| | Water content | 0 – 90 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,5177 nD | ± 0,0003 nD | 0,0001 nD |

Scope of application: Salt

The following models are particularly suitable to determine the concentration of NaCl (salt) in water and seawater. This is often used for the preparation and for the cooking of sauces, bases for pastries, the production of brines (e.g. for white cheese) and the preparation of seafood and marinades for meat. Alternatively the display can be switched to show Brix or the refractive index.

The main scope of applications is:

- Food industry
- Restaurants, and large-scale catering establishment, canteens
- Fisch farm



| Model | Scales | Measuring range | Accuracy | Division |
|----------------|---------------------------|--------------------|-------------|-----------|
| KERN | | | | |
| ORM 1NA | Salt content (NaCl) % | 0 – 28 % | ± 0,2 % | 0,1 % |
| | Salt content (NaCl) ‰ | 0 – 280 ‰ | ± 2 ‰ | 1 ‰ |
| | Spec. Gravity | 1,000 – 1,220 | ± 0,002 | 0,001 |
| | Brix | 0 – 28 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4100 nD | ± 0,0003 nD | 0,0001 nD |
| ORM 1SW | Salt content seawater | 0 – 100 ‰ | ± 2 ‰ | 1 ‰ |
| | Chlorine content seawater | 0 – 57 % | ± 2 ‰ | 1 ‰ |
| | Spec. Gravity | 1,000 – 1,070 | ± 0,002 | 0,1 % |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |

Digital refractometer KERN ORM-AL · ORM-BR · ORM-WN

Scope of application: Beer/alcohol

The following models are particularly suitable for determining the sugar content of the original wort of beer in its unfermented state. The value can be read straightaway, without having to be converted, using the SG Wort and Degrees Plato scales. In addition, the percent by volume and percent by mass scales can be used to determine the alcohol content of clear spirits.

The main scope of applications is:

- Beer brewers
- Alcohol production



| Modell | Skalen | Messbereich | Genauigkeit | Teilung |
|----------------|----------------------|--------------------|-------------|-----------|
| KERN | | | | |
| ORM 1AL | Percentage by mass | 0 – 72 % | ± 1 % | 1 % |
| | Percentage by volume | 0 – 80 % | ± 1 % | 1 % |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |
| ORM 1BR | Plato | 0 – 31 °P | ± 0,3 °P | 0,1 |
| | SG Wort | 1,000 – 1,130 | ± 0,002 | 0,1 |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |

Scope of application: Wine

The following models are particularly suitable for the measurement of the sugar content in fruit. It indicates the expected °Alcohol of the fruit. The degree of ripeness of fruit (fruit-sugar) can also be determined, such as e.g. grapes. Alternatively the display can be switched to show Brix or the refractive index.

The main scope of applications is:

- Agriculture: Wine-growing (viticulture) and fruit-growing
- Wine-production
- Must and alcohol production



°Oe = Degree Oechsle, °KMW = Klosterneuburger Most Waage

| Model | Scales | Measuring range | Accuracy | Division |
|----------------|----------------------|-----------------|------------|----------|
| KERN | | | | |
| ORM 1WN | Oechsle | 0 – 150 °Oe | ± 2 °Oe | 1 °Oe |
| | Percentage by volume | 0 – 22 % | ± 0,2 % | 0,1 % |
| | KMW (Babo) | 0 – 25 °KMW | ± 0,2 °KMW | 0,1 °KMW |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| ORM 2WN | Oechsle France | 0 – 230 °Oe | ± 2 °Oe | 1 °Oe |
| | Percentage by volume | 0 – 22 % | ± 0,2 % | 0,1 % |
| | KMW (Babo) | 0 – 25 °KMW | ± 0,2 °KMW | 0,1 °KMW |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |

Digital refractometer KERN ORM-CO · ORM-UN

Scope of application: Coffee

The following models are particularly suitable for measuring the dissolved solids (TDS) in coffee to determine or compare the strength of a cup of coffee. For roasting plants, the TDS% value is used to determine the solubility level of a roast and to control the quality. Alternatively the display can be switched to show Brix or the refractive index.



The main scope of applications is:

- Coffee industry
- Coffee roasting plants
- Coffee competitions

| Modell | Skalen | Messbereich | Genauigkeit | Teilung |
|----------------|------------------|--------------------|-------------|-----------|
| KERN | | | | |
| ORM 1CO | Coffee TDS 1 | 0 – 25 | ± 0,2 | 0,1 |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |
| ORM 2CO | Coffee TDS 2 | 0 – 25 | ± 0,2 | 0,01 |
| | Brix | 0 – 30 | ± 0,2 | 0,1 |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |

Scope of application: Urine

The following models are particularly suitable for the measurement of the specific gravity (sg) in urine, the quantity of serum (serumproteine) in urine (doping control among athletes), and the refractive index.



The main scope of applications is:

- Hospitals
- Doctor's surgeries/Physicians
- Medical training institutions
- Nursing homes
- Sports medicine (doping test)
- Veterinary

| Model | Scales | Measuring range | Accuracy | Division |
|----------------|-----------------------|--------------------|----------------|--------------|
| KERN | | | | |
| ORM 1UN | Urine (spec. gravity) | 1,000 – 1,050 sgU | ± 0,001 sgU | 0,001 sgU |
| | Serum protein | 0 – 12 g/100 ml | ± 0,2 g/100 ml | 0,1 g/100 ml |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |
| ORM 2UN | Urine (s. g. dog) | 1,000 – 1,060 sgU | ± 0,002 sgU | 0,001 sgU |
| | Urine (s. g. cat) | 1,000 – 1,060 sgU | ± 0,002 sgU | 0,001 sgU |
| | Brix | 0 – 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 – 1,4200 nD | ± 0,0003 nD | 0,0001 nD |

Digital refractometer KERN ORM-CA

Scope of application: Industry/Automotive
































The following models are particularly suitable for the measurement and determination of AdBlue®, glycol concentration (ethylene (EG) and propylene (PG)), battery fluid (BF), urea, the freezing point of windscreen wash water (CW). Furthermore these models are suitable for the measurement of thermal exchange systems. Alternatively the display can be switched to show Brix.

The main scope of applications is:

- Automotive industry: Car-workshops and producers
- Chemical industry
- Solar industry: Antifreeze monitoring



| Model | Scales | Measuring range | Accuracy | Division |
|----------------|-----------------------|-----------------|-------------|-------------|
| KERN | | | | |
| ORM 1CA | Wash water | (-60) - 0 °C | ± 0,5 °C | 0,1 °C |
| | AdBlue® | 0 - 51 % | ± 0,2 % | 0,1 % |
| | Battery fluid | 1,000 - 1,500 | ± 0,005 | 0,1 % |
| | Brix | 0 - 50 % | ± 0,2 % | 0,1 % |
| | Refractive index | 1,3330 - 1,4200 | ± 0,0003 nD | ± 0,0001 nD |
| ORM 2CA | Ethylene glycol (%) | 0 - 100 % | ± 0,5 % | 0,1 % |
| | Ethylene glycol (°C) | (-50) - 0 °C | ± 0,5 °C | 0,1 °C |
| | Propylene glycol (%) | 0 - 100 % | ± 0,5 % | 0,1 % |
| | Propylene glycol (°C) | (-60) - 0 °C | ± 0,5 °C | 0,1 °C |
| | Brix | 0 - 90 % | ± 0,2 % | 0,1 % |

- 
360° rotatable microscope head
- 
Monocular Microscope
 For the inspection with one eye
- 
Binocular Microscope
 For the inspection with both eyes
- 
Trinocular Microscope
 For the inspection with both eyes and the additional option for the connection of a camera
- 
Abbe Condenser
 With high numerical aperture for the concentration and the focusing of light
- 
Halogen illumination
 For pictures bright and rich in contrast
- 
LED illumination
 Cold, energy-saving and especially long-life illumination
- 
Incident illumination
 For non-transparent objects
- 
Transmitting illumination
 For transparent objects
- 
Fluorescence illumination
 For stereomicroscopes
- 
Fluorescence illumination for compound microscopes
 With 100 W mercury lamp and filter
- 
Fluorescence illumination for compound microscopes
 With 3 W LED illumination and filter
- 
Phase contrast unit
 For a higher contrast
- 
Darkfield condenser/unit
 For a higher contrast due to indirect illumination
- 
Polarising unit
 To polarise the light
- 
Infinity system
 Infinity corrected optical system
- 
Zoom magnification
 For stereomicroscopes
- 
Auto-focus
 For automatic control of the focus level
- 
Parallel optical system
 For stereomicroscopes, enables fatigue-proof working
- 
Integrated scale
 In the eyepiece
- 
SD card
 For data storage
- 
USB 2.0 digital camera
 For direct transmitting of the picture to a PC
- 
USB 3.0 digital camera
 For direct transmitting of the picture to a PC
- 
WiFi data interface:
 For transmitting of the picture to a mobile display device
- 
HDMI digital camera
 For direct transmitting of the picture to a display device
- 
PC software
 To transfer the measurements from the device to a PC.
- 
Automatic temperature compensation
 For measurements between 10 °C and 30 °C
- 
Protection against dust and water splashes IPxx:
 The type of protection is shown in the pictogram cf. DIN EN 60529:2000-09, IEC 60529:1989+A1:1999+A2:2013
- 
Battery operation
 Ready for battery operation. The battery type is specified for each device.
- 
Battery operation rechargeable
 Prepared for a rechargeable battery operation
- 
Plug-in power supply
 230V/50Hz in standard version for EU. On request GB, AUS or USA version.
- 
Integrated power supply unit
 Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.
- 
Package shipment
 The time required to manufacture the product internally is shown in days in the pictogram.

ABBREVIATIONS

- C-Mount** Adapter for the connection of a camera to a trinocular microscope
- FPS** Frames per second
- H(S)WF** High (Super) Wide Field (Eyepiece with high eye point for wearers of glasses)
- LWD** Long Working Distance
- N.A.** Numerical Aperture
- SLR camera** Single-Lens Reflex camera
- SWF** Super Wide Field (Field number at least \varnothing 23 mm for 10 \times eyepiece)
- W.D.** Working Distance
- WF** Wide Field (Field number up to \varnothing 22 mm for 10 \times eyepiece)



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.