Serie B-510

MANUALE D'ISTRUZIONI

	Modello
	B-510POL
	B-510POL-I

v 1.0 2018

Indice

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v 1.0 2018

Table of contents

- 1. Warning**
 - 2. Symbols and conventions**
 - 3. Safety informations**
 - 4. Intended use**
 - 5. Overview**
 - 6. Unpacking**
 - 7. Assembling**
 - 8. Summary of transmitted light brightfield observation procedures (B-510POL / B-510POL-I)**
 - 9. Use of the microscope (transmitted light - brightfield)**
 - 10. Use of the microscope (polarized light)**
 - 11. Microphotography**
 - 12. Maintenance**
 - 13. Troubleshooting**
- Equipment disposal**

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information


Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5. Overview (B-510POL)

5. Overview (B-510POL-I) (opposite side)

6. Unpacking (B-510POL)

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the B-510POL box, the microscope parts are the following:

- | | |
|------------------------------|---------------------------|
| ① Microscope frame | ⑨ Tint plates |
| ② Eyepieces | ⑩ Tension adjustment tool |
| ③ Objectives | ⑪ Allen wrench |
| ④ Observation head | ⑫ Dummy slider |
| ⑤ Bertrand lens | ⑬ Power supply |
| ⑥ Analyzer | |
| ⑦ Nosepiece centering screws | |
| ⑧ Dust cover | |

6. Unpacking (B-510POL-I)

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling


Once opened the B-510POL-I box, the microscope parts are the following:

- | | |
|------------------------------|-------------------------------|
| ① Microscope frame | ⑨ Tint plates |
| ② Eyepieces | ⑩ Tension adjustment tool |
| ③ Objectives | ⑪ Allen wrench |
| ④ Observation head | ⑫ Reflected light illuminator |
| ⑤ Bertrand lens | ⑬ Reflected light polarizer |
| ⑥ Analyzer | ⑭ Power supply |
| ⑦ Nosepiece centering screws | |
| ⑧ Dust cover | |

Assembling procedure (B-510POL)

1. Insert the Bertrand lens ① in the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 1)

2. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig.2)

- ▶ **Hold the head with one hand during the locking in order to avoid that the head falls**

3. Insert both eyepieces into the tubes of the optical head. (Fig.3)

- ▶ **One of the two eyepieces is equipped with a crosshair for centering the entire optical system. It is advisable to insert the eyepiece with crosshair in the right eyepiece holder.**

4. The condenser is pre-installed before leaving the factory. To remove the condenser use an Allen wrench diam. 1,5 and operate on the locking screw placed on the right side of the condenser holder.

5. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig.4)

Assembling procedure (B-510POL)

6. Remove the dummy slider from Bertrand lens ③ and insert the analyzer ④. (Fig. 5-6)


7. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig.7)


DO NOT DISASSEMBLE THE INSTRUMENT
Do not disassemble the instrument.
This will void the warranty and can cause malfunctions.

Assembling procedure (B-510POL-I)

1. Insert the reflected light illuminator ① in the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 8)

2. Connect the illuminator ③ to the socket ④ placed in the back side of the frame. (Fig.9)

3. Insert the Bertrand lens ⑤ in the frame and lock the locking screw ⑥ with the provided Allen wrench (Fig. 10)

4. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig.11)

- ▶ **Hold the head with one hand during the locking in order to avoid that the head falls.**

5. Insert both eyepieces into the tubes of the optical head. (Fig.12)


- ▶ **One of the two eyepieces is equipped with a crosshair for centering the entire optical system. It is advisable to insert the eyepiece with crosshair in the right eyepiece holder.**

Assembling procedure (B-510POL-I)

6. The condenser is pre-installed before leaving the factory. To remove the condenser use an Allen wrench diam. 1,5 and operate on the locking screw placed on the right side of the condenser holder.


7. Insert reflected light polarizer ⑦. (Fig. 13)

8. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig.14)

9. Remove the dummy slider from Bertrand lens ⑧ and insert the analyzer ⑨. (Fig. 15-16)


10. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig.17)

DO NOT DISASSEMBLE THE INSTRUMENT
Do not disassemble the instrument.
This will void the warranty and can cause malfunctions.

8. Summary of transmitted light brightfield observation procedures (B-510POL / B-510POL-I)

9. Use of the microscope (transmitted light - brightfield)

1. Light intensity adjustment

Operate on the light intensity adjustment knob to turn ON / OFF the microscope and to increase / decrease the illumination voltage ①. (Fig.18)

2. Coarse focus tension adjustment (Fig. 19)

► Adjust the tension using the provided tool.

The coarse knob tension is pre-set in the factory.

To modify the tension according to personal's needs, rotate the ring ② using the provided tool.

Clockwise rotation increases the tension.

If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

3. Coarse upper limit knob (Fig. 20)

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

After focussing the specimen, rotate the knob ③ and lock it. In this way the focus upper limit is set. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.

Fine focus movement is not affected by the coarse focus lock.

► To unlock, move the knob in the opposite direction to the one used for the lock.

4. Stage (Fig. 21)

The rotating stage accepts specimens on slide (B-510POL) or opaque specimens (B-510POL-I).

It is possible to lock the specimen once placed on the stage using the stage clips ①.

After loosening the locking screw ②, the stage can be rotated by 360°.

Fig.18

Fig.19

Fig.20


Fig.21

5. Dioptic adjustment (Fig. 22)

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ① to compensate. (Fig.22)

- ▶ **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptic correction.**

Fig.22

6. Adjusting the interpupillary distance (Fig. 23)

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- ▶ **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator's eyes. (Fig.23)**

The range of the interpupillary distance is 48-75mm.

Fig.23

7. Use of eye shields (Fig.24-25)

- **Use with eyeglasses**
Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses.
- **Use without eyeglasses**
Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation.

Fig.24

Fig.25

8. Centering the condenser (Fig.26)

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①.
3. Rotate the field diaphragm ring ② in the direction showed by the arrow, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circoscribes the field of view.

Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces.

Aperture diaphragm (Fig. 27)

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 28.

Example: with objective PLAN 40x / 0,65 set the scale to $0.65 \times 0.8 = 0,52$

Fig.26

CENTER THE CONDENSER


Fig.27


Fig.28

9. Centering of reflected light diaphragms

Field diaphragm (FS)

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter "R". (Fig.29)
2. Place a specimen on the stage, insert 10x objective and focus the specimen.
3. Rotate the knurled ring of the field diaphragm ① in the direction showed by the arrow to fully close the diaphragm. (Fig.30)
4. Using the provided Allen wrench use the two centering screws ② bring the bright spot in the center of the field of view.
5. Gradually open the diaphragm. The illuminator is centered when the diaphragm image is symmetrical to the field of view.
6. In normal use, open the diaphragm until it circoscribes the field of view.

Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces.

Aperture diaphragm (AS)

1. Rotate the knurled ring of the aperture diaphragm ③ in the direction showed by the arrow to fully close the diaphragm.
2. Remove one eyepiece.
3. While observing in the empty eyepiece holder, use the provided Allen wrench and use the centering screws ④ to bring the bright spot in the center of the field of view. (Fig.31)
4. The illuminator is centered when the diaphragm image is symmetrical to the field of view
 - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image
 - With low contrast specimens rotate the ring on the illuminator to set the numerical aperture value to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the ring of the illuminator in order to obtain an image like the one in fig. 28.


Fig.29

Fig.30


Fig.31




10. Use of the microscope (polarized light)

- The system allows observation in Orthoscopy (crossed Nicol) or in Conoscopy (crossed Nicol with the use of the Bertrand lens).
- For optimal performance in polarized light microscopy, accurate optical adjustments are essential before beginning the observation.

1) INITIAL SETUP

Centering of rotating stage

1. Center the condenser as described in Chapter 9.
2. Loosen the stage rotation lock ① and rotate the stage until the graduated scale of the stage ② and the vernier scale are aligned ③ in the "0" position. (Fig.32)
- **This operation serves to ensure a standard reference position for centering the rotating stage.**
3. Focus on a recognizable detail ④ in the field of view by placing it at the center of the crosshair. (Fig.33)
4. By turning the stage, the detail in focus will describe a circle ⑤. (Fig.33)
5. Move again the stage on the "0" position and lock the locking screw ①. Using the centering screws of the stage ⑥ move the detail in diametrically opposite direction to the described circle. The movement must be about half the diameter of the circle described. (Fig.34)
6. Manually move the specimen and bring it back to the center of the crosshair. Loosen the stage locking screw again and rotate the stage again.
7. If the centering has been carried out correctly, by turning the table the image of the focused detail does not move with respect to the center of the crosshair. If not, repeat the operations described by (2) to (7) until the perfect coincidence of the center of rotation of the stage with the center of the crosshair for which the preparation remains in the center of the crosshair turning the table.
8. Once centered with the 10x, rotate the nosepiece to insert in the optical path all the other objectives and verify the correct centering of the objectives, acting on the screws of centering the nosepiece ①, to make sure that all the objectives are perfectly centered with respect to the optical axis. (Fig.35)

► **B-510POL**

Checking the extinction of light (“crossed Nicol” position)

1. Remove the specimen from the light path and insert 10x.
2. Insert the condenser polarizer, loosen the locking screw of the polarizer ① and check that the scale is in the “0” position ②. (Fig.36)
3. Insert rotatable analyzer in the light path, loosen the locking screw ③ and put the vibration scale on 0° ④, then lock the locking screw ③. (Fig.37)
4. Loosen the polarizer locking screw ① and rotate the polarizer scale ② to obtain total extinction (total dark in the eyepieces). Tighten the screw ①. (Fig.36)

► **It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but is due to the mechanical alignment of the polarizers during assembly.**

► **B-510POL-I**

Checking the extinction in reflected light

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter “R”. (Fig.38)
2. Insert the reflected light polarizer ②. (Fig.39)
3. Place on the stage a flat mirror and insert 10x objective.
4. Insert rotatable analyzer in the light path, loosen the locking screw ③ and rotate the vibration direction scale ④ to obtain complete dark in the eyepieces, then lock the locking screw ③. (Fig.37)

► **It may happen that the scale of the analyzer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but is due to the mechanical alignment of the polarizers during assembly.**

Checking the extinction in transmitted light

1. Move the selector ① on the reflected light illuminator in the fully removed position, corresponding to letter “T”. (Fig.38)
2. Repeat procedure described in steps 1. to 4. for the B-510POL.


Fig.36

Fig.37


Fig.38


Fig.39

2) Use of tint plates

Three tint plates are supplied with the microscope:

-) λ plate (1st order Red)
-) $\lambda/4$ plate
-) "Quartz wedge" (Q) plate


1. Insert in the right slot of the Bertrand lens ① one of the tint plates ②. (Fig.40)
 2. In polarized light, inserting one of the plates will have chromatic effects on the specimen.
- Using the λ plate (also called 1st order Red) the specimen will take a magenta tinge.
 - Using the $\lambda/4$ the specimen will take on a color tending to pale yellow.
 - Using the Q plate the specimen will present a series of colored bands that will fade as the plate is inserted.

3) Use of Bertrand lens

Bertrand lens allows observation in Orthoscopy and Conoscopy.

In the disengaged position ("O") the lens allows observation in Orthoscopy, while in the inserted position ("B") it is possible to make observations in Conoscopy.

1. Rotate the upper knurled ring of the Bertrand lens ① to engage the "B" position. (Fig.41)
 2. Using one objective from 20x to 60x, focus the conoscopic image using the focus ring ②.
 3. If the conoscopic image is not perfectly centered with respect to the optical axis, center the image using the centering screws ③.
- By turning the stage you will see black fringes that will appear and disappear depending on the rotation of the stage. These fringes are the crystallization axes of that specific crystal.(Fig.42)

11. Microphotography

Installing the C-mount adapter

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig.43)
2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig.44)

Fig.43

Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 45).
 - "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
- **If using an SLR camera, mirror movement may cause the camera to vibrate.**
We suggest lifting the mirror, using long exposure times and a remote cord.

Fig.44

Fig.45

17. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

13. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
	Bertrand lens is in	Remove Bertrand lens from light path
	You are in a position of extinction	Disengage analyzer from the light path
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged	Make sure that the revolving nosepiece clicks properly into place.
	Tint plate or Bertrand lens are in a intermediate position.	Move to a click stop
Dirt or dust is visible in the field of view..	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	lean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings.
Visibility is poor. · Image is not clear. · Contrast is poor. · Details are indistinct. · Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
Conoscopic image cannot be seen	Swing-out lens of the condenser is not in the light path	Insert the lens in the light path
	Bertrand lens is not in the light path.	Insert the lens in the light path.
Total extinction cannot be obtained	Analyzer is not in the light path	Insert analyzer in the light path.

II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

OPTIKA® S.r.l.

Via Rigla, 30 - 24010 Ponteranica (BG) - ITALIA Tel.: +39 035.571.392 - Fax: +39 035.571.435
info@optikamicroscopes.com - www.optikamicroscopes.com

OPTIKA® Spain

spain@optikamicroscopes.com

OPTIKA® USA

usa@optikamicroscopes.com

OPTIKA® China

china@optikamicroscopes.com

OPTIKA® Hungary

hungary@optikamicroscopes.com

OPTIKA® India

india@optikamicroscopes.com