Serie B-510

MANUALE D'ISTRUZIONI

Modello
B-510BF
B-510ERGO
B-510ASB
B-510PH
B-510FL
B-510LD1
B-510LD2
B-510-2
B-510-2F
B-510-3
B-510-5

v 1.3 2018

Indice

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v 1.3 2018

Table of Contents

1. Warning
 2. Symbols and conventions
 3. Safety Information
 4. Intended use
 5. Overview
 6. Unpacking
 7. Assembling
 8. Summary of brightfield observation procedures (B-510BF/B-510ERGO)
 9. Use of the microscope
 10. Use of universal condenser for brightfield / darkfield / phase contrast (B-510PH)
 11. Microphotography
 12. Overview (B-510FL, B-510LD1/LD2)
 13. Use of the microscope in fluorescence (B-510FL, B-510LD1/LD2)
 14. Summary of fluorescence observation procedures (B-510FL)
 15. Summary of fluorescence observation procedures (B-510BLD1/LD2)
 16. Simultaneous observation in Phase contrast + Fluorescence (B-510FL only)
 17. Maintenance
 18. Troubleshooting
- Equipment disposal

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information


Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.


4. Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.


5. Overview (B-510BF/B-510ERGO)


5. Overview (B-510BF) (opposite side)

5. Overview (B-510PH)

5. Overview (B-510PH) (opposite side)

6. Unpacking (B-510BF)

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the B-510BF/B-510ERGO box, the microscope parts are the following:

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Immersion oil

- ⑥ Allen wrench
- ⑦ Tension adjustment tool
- ⑧ Dust cover
- ⑨ Power supply

6. Unpacking (B-510PH)

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the B-510PH box, the microscope parts are the following:

① Microscoper frame

② Eyepieces

③ Objectives

④ Observation head

⑤ Centering telescope

⑥ Immersion oil

⑦ Allen wrench

⑧ Tension adjustment tool

⑨ Dust cover

⑩ Green filter

⑪ Power supply


Assembling procedure

1. Insert the optical head above the stand and tighten the screw. (Fig.1)
▶ **Hold the head with one hand during the locking in order to avoid that the head falls.**
2. Insert both eyepieces into the tubes of the optical head. (Fig.2)
3. Condenser is pre-installed in the factory. To remove the condenser use an Allen wrench 1,5mm diam and operate on the locking screw placed on the right side of the condenser holder.
4. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig.3)
5. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig.4)

DO NOT DISASSEMBLE THE INSTRUMENT
Do not disassemble the instrument.
This will void the warranty and can cause malfunctions.

8. Summary of brightfield observation procedures (B-510BF)

9. Use of the microscope

1. Light intensity adjustment

Operate on the light intensity adjustment knob to turn ON / OFF the microscope and to increase / decrease the illumination voltage ①. (Fig.5)

2. Coarse focus tension adjustment

▶ Adjust the tension using the provided tool.

The coarse knob tension is pre-set in the factory.

To modify the tension according to personal's needs, rotate the ring ② using the provided tool (Fig. 6).

Clockwise rotation increases the tension.

If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

3. Coarse upper limit knob (Fig. 7)

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

After focussing the specimen, rotate the knob ③ and lock it. In this way the focus upper limit is set. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.

Fine focus movement is not affected by the coarse focus lock.


▶ To unlock, move the knob in the opposite direction to the one used for the lock.

4. Stage (Fig. 8)

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverside 0,17mm. It is possible to place two slides side by side on the stage.

- Open the spring arm of the slide holder ④ and place frontally the slides on the stage.
- Gently release the spring arm of the slide holder.
- ▶ A sudden release of the the spring arm could cause the falling of the slide.

5. Dioptric adjustment (Fig. 9)

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring ① to compensate. (Fig.9)

- ▶ **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction.**

Fig.9

6. Adjusting the interpupillary distance (Fig. 10)

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- ▶ **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator's eyes. (Fig.10)**

The range of the interpupillary distance is 48-75mm.

Fig.10

7. Use of eye shields (Fig.11-12)

- **Use with eyeglasses**
Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses.
- **Use without eyeglasses**
Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation.

Fig.11

Fig.12

8. Centering the condenser (Fig.13)

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①.
3. Rotate the field diaphragm ring ② in the direction showed by the arrow, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circoscribes the field of view.

Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces.



Aperture diaphragm (Fig. 14)

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 15.

Example: with objective PLAN 40x / 0,65 set the scale to $0.65 \times 0.8 = 0,52$

CENTERING THE CONDENSER

9. Use of oil immersion objective (Fig.16)




1. Focus the specimen with a low power objective.
 2. Lower the stage (remembering to lock the coarse upper limit knob).
- ▶ **Put a drop of oil (provided) on the area of the specimen to be observed. Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
- To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
 - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
3. Insert immersion objective.
 4. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
 5. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- ▶ **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

10. Use of universal condenser for brightfield/darkfield/phase contrast

Universal condenser provided with B-510PH allows observation in brightfield, darkfield and phase contrast.

OBSERVATION MODE	Condenser Turret position
Brightfield	BF (Fig. 17)
Darkfield	DF (Fig. 18)
Phase contrast (10x)	10/20 (Fig.19)
Phase contrast (20x)	10/20 (Fig.19)
Phase contrast (40x)	40 (Fig. 20)
Phase contrast (100x)	100 (Fig. 21)

1) BRIGHTFIELD OBSERVATION (BF)

Rotate the condenser turret to insert the “BF”position. Now repeat the steps described in the procedure “Summary of brightfield observation procedures” at pag.47.

2) DARKFIELD OBSERVATION (DF)




1. Rotate the condenser turret to insert the “DF”position.
 2. Open the aperture diaphragm.
 3. Place a specimen on the stage and focus.
 4. Observing into eyepieces raise or lower the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect.
- **Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.**
 - **“Dry” darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0,7.**
 - **Observing in darkfield, it may be necessary to raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.**

3) PHASE CONTRAST OBSERVATION (PH)

1. Center the condenser as already described at pag. 50.
 2. Rotate the condenser turret to insert the "10/20" position.
 3. Insert 10x objective into the light path.
 4. Open aperture diaphragm.
 5. Place a specimen on the stage and focus.
 6. Remove one eyepiece and insert the centering telescope. (Fig.22)
 7. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 23)
 8. Using centering screws on the condenser ①, center the phase rings to make the bright ring ② be concentric to the dark ring ③.
 9. Insert 20x objective (do not rotate the condenser turret) and check the centering of the two rings. (Fig. 25)
 10. Repeat the same operation with other objectives to check the ring centering: 40x objective – turret position "40", 100x objective – turret position "100".
 11. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
- **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
 - **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

Use of the green filter (Fig. 26)

- The green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope (Fig. 26) and begin the observation.
- For observation in brightfield or darkfield it is advisable to remove the filter from the optical path.

11. Microphotography

Installing the C-mount adapter

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig.27)
2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig.28)

Fig.27

Use of Reflex camera

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 29).
- "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
- **If using an SLR camera, mirror movement may cause the camera to vibrate.**
We suggest lifting the mirror, using long exposure times and a remote cord.

Fig.28

Fig.29

12. Overview (B-510FL)

The main microscope commands remain unchanged: only the fluorescence parts are highlighted.

12. Overview (B-510LD1/2)

The main microscope commands remain unchanged: only the fluorescence parts are highlighted.

13. Use of the microscope in fluorescence (B-510FL, B-510LD1/LD2)

This section refers exclusively to the use of the reflected light fluorescence microscope.
For transmitted light operations, refer to this manual in sections 8-9-10 from page 11 to page 17.

13.1. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

13.2. Assembling (B-510FL)

Once opened the B-510FL box, the microscope parts are the following:

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Lamp housing
- ⑥ Fluorescence illuminator
- ⑦ HBO bulbs

- ⑧ Allen wrenches
- ⑨ Tension adjustment tool
- ⑩ Dust cover
- ⑪ UV protective shield
- ⑫ Power supply
- ⑬ Power cord
- ⑭ Fluorescence power supply

13.2. Assembling (B-510LD1/LD2)

Once opened the B-510LD1/LD2 box, the microscope parts are the following:

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Fluorescence illuminator

- ⑥ Allen wrenches
- ⑦ Tension adjustment too
- ⑧ Immersion oil
- ⑨ Dust cover
- ⑩ Power supply

1. Assembling procedure

► B-510FL ONLY

1. Using the provided Allen wrenches, remove the lamp housing from the illuminator using the tightening screws ①. (Fig.30)

2. Insert the lamp housing extension tube and tighten the screws ②. (Fig. 31)


3. Reassemble the lamp housing and tighten the screws ①. (Fig. 32)

► FOR ALL MODELS

4. Insert the round dovetail socket of the illuminator ③ into the hole in the microscope body and tighten the locking screw ④. (Fig 33)

► **B-510FL ONLY**

- Disconnect all electrical cables before installing or replacing the bulb.
- The bulb has an anode and a cathode of different sizes. Respect the polarity during assembly, respecting the bulb dimensions.
- Do not touch the bulb of the lamp with bare hands to leave no traces of grease on the bulb. If this happens, clean the bulb with a soft cloth before turning on the lamp.
- The lamp has an average life of about 200-250 hours: a time counter and a voltage indicator are shown on the lamp power supply. Replace the lamp when the hour count exceeds 250 or if the voltage drops below 4.5A.


- During use, the lamp, the lamp housing and the surrounding environment become hot.
- Before replacing the lamp, switch off the power supply, disconnect all cables and wait for the lamp and the lamp housing to cool.
- After switching on the lamp, wait at least 10-15 minutes before switching it off.
- After switching off the lamp, wait for 5-10 minutes before switching it on again so that the mercury vapors have time to condense.
- The lamp contains ultraviolet radiation that could be harmful to eyes and skin. Always look at the lamp arc through the provided orange screen.

► **B-510FL ONLY**

1. Open the lamp housing using the door lock screw ① and remove the lamp holder. (Fig.34)

2. Remove the plastic block ② from the lamp holder (or the exhausted lamp in case of replacement) by loosening the two locking screws ③. (Fig.35)


3. Insert the mercury bulb ④ (respect the polarity of the bulb), tighten the locking screws and refit the lamp holder inside the lamp housing. (Fig. 36)

4. Insert the lamp housing cable into the fluorescent power supply, aligning the notches on the connectors. (Fig.37)


5. Insert the power cord into the connector ①. (Fig.38)

► **B-510LD1/LD2 ONLY**

6. Connect the LED illuminator plug ① into the microscope body. (Fig.39)





2. Setting the microscope

► B-510FL ONLY


Centering the HBO mercury bulb.

► Wait around 5 minutes before proceeding with this operation to allow the bulb to warm up properly.


1. Turn on the mercury vapor bulb by operating the power supply switch ①. (Fig.40)
2. Turn the nosepiece into an empty position (without objectives) and remove the protective cap, or remove an objective from the nosepiece.
3. Place a piece of white paper on the table and insert the fluorescent cube "B" into the optical path. (Fig.41)
4. Acting on the focus screw of the collector lens ① and on the centering screws ② try to obtain the light spot of the bulb's arc. (Fig.42-43)


- Using the focus screw of the collector lens ①, put the image of the arc projected onto the paper. The light spot must be brighter and sharper as possible. (Fig. 44)

- Using the centering screws ② on the side of the lamp housing, center the image of the arc. (Fig.44-45)

- Using the focusing screw of the collector lens ① enlarge the image until a homogeneous illumination is achieved. (Fig.46). At this point, insert an objective into the optical path and, looking into the eyepieces, optimize the illumination always using the screws ① and ②.

- After replacing the exhausted bulb, reset the time counter on the power supply by pressing the "Reset" button ①. (Fig.47)

3. Use of the microscope

► B-510FL ONLY

1. Turn on the power supply for the mercury bulb and wait 5 minutes for the arc to stabilize. (Fig.48).
2. Move the filter selector ④ to one of the 4 available positions until the click stop. (Fig.49)
3. The microscope has a 4-position filter holder slider. The positions 1 and 2 are empty to house additional filters, position 3 houses a B filter and position 4 a G filter.

Fig.48

FILTER CUBE	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
B	475/30 nm	505 nm	515LP nm	<ul style="list-style-type: none"> • FITC: fluorescent antibodies • Achridine orange: DNA, RNA • Auramine
G	530/40 nm	570 nm	590LP nm	<ul style="list-style-type: none"> • Rhodamine, TRITC: fluorescent antibodies • Propidium iodide: DNA, RNA • RFP

► B-510LD1/LD2 ONLY:

1. Turn on the fluorescence LED, using the switch on the back side of the microscope.
2. Move the filter selector ④ to one of the available positions until the click stop. (Fig.49).
3. LD1 and LD2 models have a 2-position filter holder slide. In the case of the LD1 model, the slide only houses a B filter, while in the LD2 model the slide houses a B and a G filter.

Fig.49


FILTER CUBE	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
B	475/30 nm	505 nm	515LP nm	<ul style="list-style-type: none"> • FITC: fluorescent antibodies • Achridine orange: DNA, RNA • Auramine
G	530/40 nm	570 nm	590LP nm	<ul style="list-style-type: none"> • Rhodamine, TRITC: fluorescent antibodies • Propidium iodide: DNA, RNA • RFP

Use of the shutter:


- ▶ **The microscope is equipped with a shutter ① located on the right side of the fluorescent illuminator. (Fig.50)**
1. Close the shutter by interrupting the observation for a limited time and not subjecting the sample to unnecessary lighting in the period in which it is not observed.
(Switching off and switching on frequently the HBO lamp considerably reduces its duration).
- ▶ **This precaution is not necessary in the case of the LD1 and LD2 models: the LED can be switched on and off without any problem.**

14. Summary of fluorescence observation procedures (B-510FL)

15. Summary of fluorescence observation procedures (B-510LD1/LD2)

16. Simultaneous observation in Phase Contrast + Fluorescence (B-510FL only)

- **This microscope allows observation in transmitted light Phase contrast in combination with reflected light Fluorescence. Samples with rapid decay must first be observed in Fluorescence and then in Phase Contrast. The combined observation allows you to easily identify some areas of the sample that emit fluorescence.**
1. Turn on the power supply for the HBO fluorescent bulb and wait 5 minutes before the arc stabilizes.
 2. Move the filter selector to an empty position or, if the filter holder is complete, to the position containing the UV filter.
 3. Insert the desired PH lens and rotate the phase contrast condenser turret to the position containing the corresponding phase ring.
 4. Focus the sample.
 5. Adjust the light intensity of the transmitted light.
 6. Move the fluorescence filter selector to the desired position.
 7. To obtain the proper observation of the sample, adjust the light intensity of the transmitted light to modulate the intensity of the fluorescence with that of the phase contrast.

17. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

18. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Fluorescence filter selector is not in a click stop	Move the selector to a click stop
	Fluorescence shutter is closed	Open the shutter to shutter
	Fluorescence filter is not suitable for the specimen	Use a suitable filter
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	The turret of the phase contrast condenser is in an incorrect position	Move the turret to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings.
Visibility is poor. · Image is not poor. · Contrast is poor. · Details are indistinct. · Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	For phase contrast observation, a brightfield objective is used instead a phase contrast one	Use a phase contrast objective
	Phase rings of objective and condenser are not well centered	Operate on centering screws
	Objective in use is not compatible with condenser phase ring	Use a compatible objective
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality

II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.