 Operating Instructions

 Mode d'emploi

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 Bruksanvisning

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 Bruksanvisning

 Bruksanvisning

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Introduction

Thank you for choosing the Prestige Medical 2100 Classic Autoclave designed to sterilize solid unwrapped instruments. After removing from the box, please check for any transit damage. If any damage to the box is found, please contact your supplier immediately

Together with this unit and operating manual, you will find the following:

- Electrical mains cord
- Instrument furniture
- Spare fuses x2 (230v models only)
- Performance test certificate
- Certificate of Compliance (230v models only)
- Warranty registration card

Key to Pictures, Diagrams and Symbols

The following descriptions refer to the pictures of the controls, display lights and operating symbols on Page 3 of this manual.

Controls:

- A - Depressurization Valve
- B - Start Cycle Button
- C - Pressure Rise Indicator
- D - Display Panel
- E - Pressure Gauge ("Plus" models only)
- F - Temperature Gauge ("Plus" models only)
- G -Top Handle - WARNING HOTPART! –
Do Not Touch the top Handle immediately after a cycle
Do not use to open the lid.
Do not carry the autoclave when HOT.
- H -Steam Duct – WARNING HOTPART!
Do Not Cover

Figures:

- F1 -Lid aligned with body, turn clockwise to close
- F2 -Lid in closed position

Warning Symbols:

- W1 - Caution, electric shock hazard
- W2 - Warning, read manual before using the autoclave
- W3 - Warning, unit must be earthed/grounded
- W4 - Warning, heat hazard

Display Lights:

- L1 -Power "On" light - illuminates GREEN
- L2 -Heating Light - illuminates ORANGE
- L3 -Sterilizing light - illuminates YELLOW
- L4 -Sterilizing complete light - illuminates GREEN
- L5 - Fault light - illuminates RED

WARNING!

Do not touch the body or lid as these parts become hot when the autoclave is in operation

Getting Started

Before using the autoclave for the first time please take time to read the following pages to familiarise yourself with the operation of the unit. We strongly recommend that all users of the autoclave are trained in its operation. The autoclave is very easy to use. By following this simple operating sequence in conjunction with the pictures of the autoclave, it's controls, display panel and operating symbols (page 3), you will be able to ensure your instruments are correctly sterilized every time.

WARNING!

Mains outlet MUST BE EARTHED (GROUNDED)

The mains plug should always be easily accessible as it is to be relied upon as "the means of disconnection"

1. Water

Fill the unit to the water level line on the inside of the chamber with 0.75l of distilled or de-ionised water.

DO NOT USE TAP WATER OR OVERFILL.

2. Loading

Place UNWRAPPED and WASHED instruments ONLY plus a TST Strip into the instrument basket or cassettes, and place in the unit.

TST strips should be placed as near to the centre of the load as possible (Available as an accessory).

-Sterilizing Cycle is completed when:

The Buzzer sounds. Lights:

L1 illuminates GREEN

L4 illuminates GREEN

Note: L4 remains GREEN until a new cycle is started or the unit is disconnected from the mains power.

6. Depressurizing.

Once the Sterilizing Cycle has been completed, the unit needs to be depressurized and allowed to cool down before the lid and sterilized instruments can be removed. The time taken to reach the point at which this is safe to do so can be shortened by Manually Depressurizing the unit.

Open the Depressurization Valve (A) by turning slowly in an anti clockwise direction*. The Pressure Indicator (C) will drop once the steam has been released.

* Warning: There will be a visual and audible release of steam from the rear of the Steam Duct (H).

7. Unlocking.

Once the pressure has been released the lid may be unlocked. Ensure the Depressurization Valve (A) is open. Remove the lid by turning in an anti-clockwise direction (F1).

8. Unloading.

Lift off the lid, gently place upside down on a solid work surface and leave to cool. Ensure the Depressurization Valve (A) is in the closed position to avoid damaging it.

The unit has completed a successful cycle if the "spot" on the TST Indicator Strip has completely changed colour from yellow to purple.*

The basket or cassettes containing the sterilized instruments can now be lifted out of the unit.

To avoid damage, replace the lid as described in step "3"

* Please Note: If the "spot" has not completely changed colour, replace with a new TST strip and start a new cycle. If the "spot" fails to change colour for a second time, **do not use the unit** until it has been checked by a qualified engineer.

STERILIZATION CYCLE HAS NOT BEEN ACHIEVED.

WARNING!

The manufacturers of instruments should be consulted about their suitability for autoclaving and the maximum temperature which the instruments can withstand.

3. Closing.

Always place the lid on the autoclave with the Depressurization Valve (A) open. Align the lid and body (F1) and turn in a clockwise direction ensuring the lid is completely closed (F2). Close the Depressurization Valve (A) so it is aligned with the top handle on the lid. Never leave in position as shown in F1.

4. Power Connection.

Attach the cable supplied to the rear of the unit and plug into an EARTHED mains electrical socket of the CORRECT voltage. Lights: L1 illuminates GREEN

5. Starting.

Start the sterilizing cycle by pressing button (B).

Lights: L1 illuminates GREEN
L2 illuminates ORANGE

-as the temperature rises, air will be displaced by steam through the Air Bleed Device located in the lid, until it closes with an audible "click", sealing the unit.

The Pressure Indicator (C) will rise indicating the unit is now pressurised.

-Sterilizing Temperature is reached when:

Lights: L1 illuminates GREEN
L2 flashes ORANGE
L3 illuminates YELLOW

9. Biological Monitoring

9.1 Selecting Biological Indicator

Healthcare personnel should select Biological Indicators consisting of spores of *Bacillus Stearothermophilus* that comply with the American National Standard, Biological Indicators for saturated steam sterilization processes in healthcare facilities (AAMI, 1986). In addition, data should be obtained from manufacturers on the reliability, safety and performance characteristics of their products. Manufacturers of Biological Indicators should also be required to provide written instructions on the storage, handling, use and microbiological testing of their products.

9.2 Frequency of Use of Biological Indicators

Table top sterilizers should be biologically monitored during installation and after any major repairs. In addition, sterilization loads should be biologically monitored at least once a week, but preferably daily. Each load containing implantable devices should be monitored and, whenever possible, the implantable devices should be quarantined until the results of the Biological Indicator testing are available. Biological Indicator should also be used for periodic monitoring of all types of packages and trays processed.

9.3 Biological Indicators

This sterilizer is designed for surface sterilization, therefore, a spore strip or self-contained Biological Indicator offers an appropriate challenge to measure delivery of an adequate dose to kill spores. Placement of Biological Indicators should be near the centre of the instrument load in the basket or in one of the cassettes. For small loads, the indicator should be placed amongst the instruments which are being processed.

9.4 Routine Biological Monitoring

For routine biological monitoring, the sterility of the load is evidenced by the killing (failure to recover) of all spores on the test Biological Indicators (spore strips). All Biological Indicator results, including results from controls, must be interpreted by a qualified individual and must be included in the sterilizer records.

Continued Operation

To ensure your autoclave gives you the years of service for which it was designed, it is important to remember a few “do’s” and “don’t’s” with regards to the operation of the unit and to carry out simple care and maintenance procedures on a weekly basis.

Do ensure that....

- 1....you read these instructions and always follow the operating sequence.
- 2....the work surface on which you will place the autoclave is flat, solid and heat resistant.
- 3....the instruments are designed to withstand the sterilizing temperatures selected, are thoroughly cleaned and rinsed before sterilizing, and are not any longer than the length, or exceed the load weight, specified - see “Technical data” section.
- 4....the water level is maintained regularly with clean distilled or deionised water only.
- 5....the unit is in a “draught free” environment and is positioned not less than 250mm from adjacent walls.
- 6....you only use green sealing gasket (219500) and that it is changed at the end of it’s life, if visibly damaged, or when shrinkage has occurred, see “Fault mode - 5”.
- 7....the lid is securely closed when the unit is not in use, to avoid the risk of accidental damage. Never leave in position as shown in F1.
- 8....you quote your model details, serial number and date of purchase when contacting Prestige Medical or your supplier.

Do not....

- 1....touch the unit whilst in operation - it gets HOT.
- 2....attempt to remove the lid during operation.
- 3....lose this operating instruction manual
- 4....add any chemicals whatsoever to the water.
- 5....attempt to sterilize volatile substances toxic materials or inappropriate loads.
- 6....place the unit on heat sensitive surfaces polished wood or glass.
- 7....open the Depressurization Valve (A) during the sterilization cycle.
- 8... leave the Depressurization Valve (A) in the “open” position when placing the lid upside down on a work surface.
- 9... immerse the unit or electrical cord in water when cleaning.
10. use abrasive materials or lubricants when cleaning.
- 11...drop or abuse the unit.
- 12...use in areas of risk associated with flammable materials or gasses.
- 13...attempt to change fuses until the unit has been unplugged from the mains. Only qualified persons should change fuses.
- 14...reach over the unit when removing cover, to do so may cause burns from rising heat and steam.
- 15...press start button once cycle has been started as this will re-set the cycle timer to zero.

Care and Maintenance.

WARNING!
Disconnect the autoclave from the mains power supply before cleaning.

Green Sealing Gasket.

1. Remove from inside the lid and clean with warm, soapy water.
2. Rinse thoroughly, shake dry, do not wipe.
3. Replace in the lid by tucking evenly under all lugs starting at the Gasket Offset Device. It may appear slightly wrinkled until used.
4. Replace gasket when it begins to show signs of leakage.

Autoclave.

- 1 If a new gasket leaks, or a persistent leak develops, gently clean the sealing surface of both the lid and body of the unit with a “Scotchbrite” scrubbing pad making sure you do not remove any metal. Rinse both surfaces but do not dry.
- 2 Clean both interior and exterior with warm, soapy water ensuring the electrical parts are kept dry.
- 3 Monitor the first cycle of the day to check the Air Bleed Device, which is located inside the lid, audibly “clicks” shut.
- 4 Prestige Medical recommends that your unit is calibrated at six monthly intervals.
- 5 Lubricate underside of body lugs with “vaseline” if the lid becomes stiff.

DO NOT LUBRICATE GASKET

Troubleshooting

In the event of a fault occurring during any stage of the units operation, identify the fault by referring to the descriptions below. The fault can be rectified by following the Fault Remedy applicable to the problem incurred.

Fault Indication/Description/Remedy

Fault 1: No Power to Unit

Light: L1 fails to illuminate.
Blown fuse / Defective Socket
/ Mains not connected

Ensure mains lead is connected.
Check / replace fuse: Check
power to socket:

Fault 2: Low water or Boil dry.

Light: L5 flashes RED

Allow unit to cool before refilling to the correct level.
Disconnect from mains then reconnect and repeat cycle.
If the fault repeats with sufficient water, arrange for a
service engineer to visit.

Fault 3: Sterilization failed to be achieved.

Light: L4 fails to illuminate GREEN and there is no audible
buzzer

Disconnect from mains then reconnect and repeat cycle.
If the fault repeats arrange for a service engineer to visit.

Fault 4: Incomplete Sterilization cycle

TST strip fails to change / completely change colour.

Check expiry date of TST strips. Disconnect from mains
then reconnect and repeat cycle. If the fault repeats
arrange for a service engineer to visit.

If the unit has been allowed to fully cool down, it is recommended that, prior to use, a warm-up cycle is run. During the course of operation the water level must be maintained up to the water level line.

Fault 5: Steam or water leaks from under the lid

i) Worn or dirty gasket.

Wash gasket and sealing surfaces on the body and lid as
described under "Care and Maintenance". If the fault
persists, replace with a new gasket.

ii) Incorrectly closed lid.

Ensure the unit is fully depressurized by opening
Depressurization Valve (A). Remove lid and re-fit carefully.
Disconnect from mains, reconnect and repeat cycle.

Fault 6: Excessive steam or water leaking from

Depressurization Valve (A).

Depressurization Valve (A) in "OPEN" position.
Close Depressurization Valve (A).

Fault 7: Incorrectly reading gauges.*

* Only applicable to "Plus" models. Disconnect from
mains then reconnect and repeat cycle. If the fault
repeats arrange for a service engineer to visit.

Technical Specification

In the unlikely event that something should go wrong, we have incorporated a number of safety features to ensure that your autoclave remains safe at all times.

Safety features.

1. Located to the rear of the lid, beneath the cover, is a spring called the Gasket Offset Device (GOD Spring), designed to prevent pressure building up if the lid has been incorrectly fitted. **DO NOT TAMPER WITH THIS SAFETY DEVICE**
2. If for any reason, the temperature falls below the minimum required sterilizing temperature, resulting in the Sterilizing Temperature light (L3) switching off, the cycle timer will re-start from zero once the correct temperature has been restored.
3. If there is an electrical or electronic failure resulting in a build up of pressure - in excess of normal operating pressure - one or all of the following safety features will be activated.
 - i) Depressurization Valve (A) will loudly and rapidly "vent" steam. On 134°C model only the overpressure Safety Valve located underneath the handle (G) will loudly and rapidly vent steam allowing its exit via the Steam Duct (H)
 - ii) The gasket will "extrude" through the slot in the rear of the lid rapidly releasing excess pressure and steam. 
 - iii) A non - resettable thermal fuse located in the base of the unit will "melt" at a pre-determined temperature, disconnecting the power.

Should any of the devices listed activate, please observe the following steps:

- a) **Do not touch the unit.**
- b) **Switch off at the wall socket and un-plug**
- c) **Allow temperature and pressure to drop before**
 - i) **touching the unit**
 - ii) **removing your instruments**
- d) **Do not attempt to re-start the unit**
- e) **Arrange for an immediate service.**

Technical Specifications - Standard Body Autoclaves & 210006

Height.....	335 mm	Capacity.....	9 litres	Max' Instrument length.....	228 mm
Width.....	340 mm	Internal chamber		Max' Load Weight.....	3.0 kilos
Net Weight.....	4.5 – 5.6 kilos	Dimensions (d/h).....	210/230 mm		

- 9 litre -134°C Body Autoclaves

Height.....	335 mm	Capacity.....	9 litres	Max' Instrument length.....	228 mm
Width.....	340 mm	Internal chamber		Max' Load Weight.....	3.0 kilos
Net Weight.....	5.6 kilos	Dimensions (d/h).....	210/230 mm		

-Extended Body Autoclaves

Height.....	420 mm	Capacity.....	12 litres	Max' Instrument length.....	290 mm
Width.....	340 mm	Internal chamber		Max' Load Weight.....	4.0 kilos
Net Weight.....	4.5 kilos	Dimensions (d/h).....	210/270 mm		

Fuses - Located under the control module, fuses F10A, 32 x 6.3mm, ceramic sand filled (230v model only), Mains plug top fuse (User replaceable), F13A to BS1362 UK ONLY.

Rating - Models are rated continuously for intermittent use.

Body - Deep drawn aluminium.

Lid – Aluminium 121/126 °C models, 304 Stainless Steel on 134°C & ASME (21006 FDA 510(K) Number K962903) units only.

Heater - Externally surface mounted mechanically fixed electric.

Temperature cut out - thermal fuse.

Pressure - calibrated pressure relief valve.

Max. Single Fault Temperature - 133.3°C (141.0°C on the 134°C model)

Over Voltage Category - Group II

Pollution Degree - Group 2

Environment Conditions - indoor use - temperature 5°C to 40°C - altitude up to 2000m - maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C. - mains supply voltage fluctuations not to exceed +10% of the nominal voltage Input Connections - Mains inlet socket 'hot' format conforming to IEC 302.

Safety Shut Down - See 'Temperature Cut Out'.

Packaging - All packaging materials are recyclable.

Additional Information

Spares

Only those spare parts supplied or specified by Prestige Medical should be used in the maintenance of the autoclave. Use of unauthorised parts will invalidate any warranty given and may adversely affect the performance and safety of the unit.

Accessories

A range of accessories are available for your autoclave as described below and pictured on page 70. Contact your supplier for full details.

- 1 - 219294 - Lifting device
- 2 - 219708 - Standard Basket
- 3 - 219291 - Cassette box
- 4 - 219706 - Cassette Rack
- 5 - 219720 - Extended Basket
- 6 - 219500 - Green Silicone Sealing Gasket
- 7 - 219277 - TST Indicator strips
- 8 - 219258 - Cord Set UK
- 9 - 219500 - Green Silicone Sealing Gasket
- 10 - 219299 - Cord Set UL

Warranty

Prestige Medical shall, in the first 12 months from the date of purchase, repair or replace free of charge any parts which prove to be defective in workmanship and/or materials. The heating element (only) is covered by a lifetime guarantee.

Prestige Medical shall not be so liable in the event that the purchaser has failed to adhere to the instructions contained herein or if the autoclave has been abused, interfered with, altered, repaired or serviced by any unauthorised party this may also result in the protection provided by the equipment being impaired.

This warranty excludes the gasket, all internal furniture and consumables. Consumer's statutory rights are not affected.

Product decontamination.

Should the unit require repair, it must be decontaminated in accordance with a recognised procedure prior to return or on-site repair. A statement of equipment contamination status must be available with the product. (Details of a suitable procedure are available on request).

Approvals.

Note: Please ask for details of approvals specific to your model.
BS3970 Part 4 (pressure vessel only).
CE Mark - Medical Device Directive (93/42/EEC)
IEC 61010 – Product safety.
BS 5500 – Pressure vessel.
ASME 8 Division 1 – Pressure vessel (120v & 134 models).
FDA 510(K) – Number K962903 (120v model).
CSA / UL – 120v model.
TSSA Canadian Design Registration – 120v model.

Clean materials

- Mild washing up liquid.
- Non abrasive cream cleaner.
- Disinfectant diluted in water – non chlorine base.

Packaging

All packaging materials used are recyclable, please dispose of accordingly.

WEEE Statement

(Waste, Electrical and Electronic Equipment)

The WEEE directive places an obligation on all EU-based manufacturers and importers to take-back electronic products at the end of their useful life.

Prestige Medical Limited accepts its responsibility to finance the cost of treatment of redundant WEEE in accordance with the specific recycling requirements.

The symbol (shown below) is present on all Prestige Medical products, which indicates that the product must NOT be disposed of with other waste. Instead it is the user's responsibility to dispose of their waste electrical and electronic equipment by handing it over to an approved reprocessor, or by returning it to Prestige Medical for reprocessing. For more information about where you can send your waste equipment for recycling, please contact your local city office or Prestige Medical.

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Prestige Medical Ltd