
Instructions

# Heat-On™ & PTFE Safety Covers

Simply the safest, fastest, most efficient way of heating a round bottom flask

Your Local Distributor

# Index

Page 1	Index
Page 2	Introduction Warranty
Page 3	Safety information Important WARNINGS
Page 4	Important WARNINGS - continued
Page 5	Products and accessories
Page 6	Products and accessories - continued
Page 7	Select your stirring hotplate Do you need an adapter plate? Optimising the temperature performance of Heat-On
Page 8	Using a separate temperature probe Optimising the stirring performance of your Heat-On Wide choice of magnetic stirring bars
Page 9	Beginning your experiment
Page 10	Multi-well formats for smaller flasks Anodised blocks or fluoropolymer coated
Page 11	Fitting optional safety lifting handles to Heat-On blocks Using a retort stand
Page 12	Flask selection process
Page 13	Set-up of PTFE Safety Covers Stirring hotplate specifications and selection guide
Page 14	WARRANTY - EMAIL BACK

**Thank you for purchasing Heat-On.**

**Please read this instruction manual thoroughly before operating your unit.**

## Introduction

Heat-On is a modular heating and stirring work station designed to accept either one or two round bottom flasks from 10ml to 5 litres. The optional PTFE Safety Covers are designed to complement and enhance the performance and safety of Heat-On blocks.

Heat-On offers the following features and benefits:

- All Heat-On blocks, inserts and the Multi-well holders, are manufactured from solid aluminium providing excellent, even heat transfer to the entire block, preventing hot spots
  - Sleek weight saving design reduces the mass of the block to be heated, improving heat-up times
- Uniquely shaped deep wells maximise heated surface area to improve heat-up time and minimise the difference between the block and solution temperature
  - Optional fluoropolymer coating gives added durability
- Tests have shown that a 250ml flask of water can be boiled in under 9 minutes, 50% faster than the competition
  - All Heat-On blocks have a probe hole for use with separate temperature probe
- Heat-On can be used up to 200°C without damage to the block. Temperatures up to 260°C can be achieved but some minor surface degradation may occur

### \*PTFE Safety Covers

These innovative, solid PTFE Safety Covers reduce the risk of users touching the 'hot' Heat-On block and provide the added benefit of lowering energy consumption.

- PTFE insulation reduces the temperature of contact surface areas by up to 50% and provides a thermal safety barrier to the operator when working at elevated temperatures
  - PTFE insulation reduces energy consumption required to heat the Heat-On block by up to 15%
- PTFE covers hot metal surfaces and helps prevent accidental spills of solvents on to hot surfaces
  - Chemically resistant PTFE withstands temperatures up to 260°C
    - Fits easily over existing Heat-On blocks
- All Heat-On Safety Covers have a probe hole for use with separate temperature probes

## Warranty

Heat-On includes one year full parts and labour warranty from date of original purchase. Warranty will only be valid if a completed **Warranty Email Back** is returned within 1 month of date of purchase (see last page).

In the event of product failure please contact your local distributor. Please do not return any goods without prior agreement.

\*This product is not available in the USA

## Safety Information

The following symbols are intended to assist the user in the safe and efficient operation of Heat-On.

	<p><b>Warning</b> Applies when there is a possibility of personal injury.</p>
	<p><b>Important Note</b> Alerts the user to important facts.</p>

## Important WARNINGS

### 1. Please read these instructions completely before using your Heat-On block.

Heat-On should only be operated by trained and competent personnel. As with all chemistries, a full risk assessment should be performed prior to starting, and care should be taken to monitor your reactions at all stages. During operation Heat-On should not be left unattended unless in a supervised area.

### 2. Operate only in a fume cupboard with protective safety sash.

### 3. During and after heating take care not to touch the Heat-On block.

### 4. Risk of burns from hot blocks

Like oil baths, Heat-On blocks remain hot for some considerable time after the heating source has been switched off. A temperature probe or temperature sensitive label can be used to indicate when blocks are too hot to touch. The use of Heat-On PTFE Safety Covers will also reduce the risk of users touching the 'hot' block.

### 5. Flask selection process

Heat-On is designed so that oversized flasks cannot fit through the slightly narrower opening at the top of the block well. This means that all suitable flasks should sit at the bottom of the block well, be a loose fit and as such spin freely in the Heat-On well.

When a flask is inserted into the well it is important to ensure that the flask touches the bottom of the well and does not 'catch or stick' on the sides of the block well. If the flask is 'catching or sticking' it may be deformed and may result in jamming and even cracking of the flask when the Heat-On block is cooling.

There is a simple rule when using glassware with Heat-On: "If it does not spin, then throw it in the bin". If the flask does not spin freely then it is unsuitable for use in the Heat-On and an alternative flask should be used.

During the design of Heat-On, Radleys evaluated and tested flasks from a wide range of glass manufacturers to ensure that Heat-On blocks will accept the widest possible range of sizes and styles. However, because the tolerances on mass produced glassware can vary greatly it is not possible to guarantee that all brands will fit. Indeed, experience shows us that glassware from a single manufacturer may vary considerably in size and shape from batch-to-batch and also from item-to-item. Customers may have to pre-select appropriate sized glassware to use with Heat-On.

### 6. Starting and stopping your reaction

In an ideal situation you should load your glassware into the block when the block is at ambient (room temperature).

However in certain situations this may not be practical. In such circumstances you should note:

- You can use cold water to cool the block down rapidly, but take care with blocks above 100°C as steam is a potential hazard.
- It is generally safe to insert flasks into blocks that are below 50°C, but the flask should be removed before the block temperature drops below the start temperature as there is a small risk that the contraction of the block could crush the glass flask.
- In all eventualities flasks should be removed at or above the temperature they were originally inserted; i.e. in the case of a flask that is inserted at 80°C, it should be removed when or before the block cools to 80°C.

### 7. Chemical resistance

Anodised Heat-On blocks are resistant to the majority of solvents and splash resistant to dilute acids and alkalis at room temperature. The fluoropolymer coated Heat-On blocks offer superior resistance to the majority of solvents, acids and alkalis.

Extended exposure to acids or alkalis will attack the surface of the Heat-On. The stronger the concentration and the longer the exposure time the more chance and intensity of any attack. Heat will also speed up and intensify any attack; e.g., hot concentrated sulphuric acid will damage the fluoropolymer coated Heat-On over a period of time. It is important to clean off any residual chemical spills immediately after they occur.


## Important WARNINGS - continued

### 8. Temperature range

Heat-On blocks can be used up to 200°C without damage to the block. Temperatures up to 260°C can be achieved but for polymer coated blocks there may be some coating degradation, while anodised blocks become discoloured.

### 9. Use of the optional PTFE Safety Covers

- The Heat-On Safety Covers are designed to be used in conjunction with Heat-On blocks. Please ensure that you have read the Heat-On instructions fully before using the Heat-On Safety Covers.
- During and after heating take care not to touch the Heat-On Safety Covers and blocks.

- Risk of burns from hot Safety Covers and blocks

Whilst the Safety Cover significantly reduces the temperature of exposed surfaces, it will still be hot and may also increase the time taken for the Heat-On block to cool down after use. A temperature probe or temperature sensitive label can be used to indicate when blocks and Safety Covers are too hot to touch.

- Chemical resistance

Heat-On Safety Covers are manufactured from solid PTFE, providing chemical resistance to the majority of chemicals. However, although the Heat-On Safety Cover is designed to give additional chemical protection to the Heat-On block, it does not provide a complete sealed enclosure, and care should be taken to avoid any spillage that may result in damage to the block.

- Temperature range

Heat-On Safety Covers can be used on Heat-On blocks operating at up to 200°C without damage to the cover or block. Operation with Heat-On block temperatures up to 260°C can be achieved but some distortion to the cover and minor surface degradation to the Heat-On block may occur.

# Products and accessories

Cat No	Description
<b>Heat-On Multi-well Holder with Fluoropolymer Coating</b>	
RR61005	Heat-On Multi-well Holder (Polymer Coated)
RR63005	PTFE Safety Cover - for Heat-On Multi-Well Holder
RR61010	Heat-On 10ml Insert (Polymer Coated)
RR61015	Heat-On 25ml Insert (Polymer Coated)
RR61020	Heat-On 50ml Insert (Polymer Coated)
RR61025	Heat-On 100ml Insert (Polymer Coated)
RR61030	Heat-On 150ml Insert (Polymer Coated)
RR61065	Heat-On Insert for 8 x 16mm Tubes (Polymer Coated)
RR61070	Heat-On Insert for 4 x 20mm Tubes (Polymer Coated)
RR61075	Heat-On Insert for 4 x 24mm Tubes (Polymer Coated)
<b>Heat-On Multi-well System with Fluoropolymer Coating</b>	
RR61000	Heat-On Multi-well Block System (Polymer Coated) Consisting of the following: 1 x RR61005 Heat-On Multi-well Holder (Polymer Coated) 2 x RR61015 Heat-On 25ml Insert (Polymer Coated) 2 x RR61020 Heat-On 50ml Insert (Polymer Coated) 2 x RR61025 Heat-On 100ml Insert (Polymer Coated) 1 x RR61085 Heat-On Adapter for 135mm Hotplates (IKA etc)
<b>Heat-On Multi-well Basic System with Fluoropolymer Coating</b>	
RR61001	Heat-On Multi-well Block System (Polymer Coated) Consisting of the following: 1 x RR61005 Heat-On Multi-well Holder (Polymer Coated) 1 x RR61015 Heat-On 25ml Insert (Polymer Coated) 1 x RR61020 Heat-On 50ml Insert (Polymer Coated) 1 x RR61025 Heat-On 100ml Insert (Polymer Coated) 1 x RR61085 Heat-On Adapter for 135mm Hotplates (IKA etc)
<b>Heat-On Blocks with Fluoropolymer Coating</b>	
RR61031	Heat-On 100ml Block with Flask Sidearm Cutouts (Polymer Coated)
RR63035	PTFE Safety Cover – for Heat-On 100ml Block
RR61032	Heat-On 200ml Florentine Block (Polymer Coated)
RR61040	Heat-On 250ml Block (Polymer Coated)
RR61041	Heat-On 250ml Block with Flask Sidearm Cutouts (Polymer Coated)
RR61042	Heat-On 300ml Florentine Block (Polymer Coated)
RR63040	PTFE Safety Cover - for Heat-On 200-300ml Block
RR61045	Heat-On 500ml Block (Polymer Coated)
RR61046	Heat-On 500ml Florentine Block (Polymer Coated)
RR63045	PTFE Safety Cover - for Heat-On 500ml Block
RR61050	Heat-On 1 Litre Block (Polymer Coated)
RR61051	Heat-On 1 Litre Florentine Block (Polymer Coated)
RR63050	PTFE Safety Cover - for Heat-On 1 Litre Block
RR61055	Heat-On 2 Litre Block (Polymer Coated)
RR61056	Heat-On 2 Litre Florentine Block (Polymer Coated)
RR61060	Heat-On 3 Litre Block (Polymer Coated)
RR61061	Heat-On 4 Litre Block (Polymer Coated)
RR61063	Heat-On 5 Litre Block (Polymer Coated)

## Products and accessories - continued

Cat No	Description
<b>Multi-well Holder with Anodised Finish</b>	
RR62005	Heat-On Multi-well Holder
RR63005	PTFE Safety Cover - for Heat-On Multi-Well Holder
RR62010	Heat-On 10ml Insert
RR62015	Heat-On 25ml Insert
RR62020	Heat-On 50ml Insert
RR62025	Heat-On 100ml Insert
RR62030	Heat-On 150ml Insert
RR62065	Heat-On Insert for 8 x 16mm Tubes
RR62070	Heat-On Insert for 4 x 20mm Tubes
RR62075	Heat-On Insert for 4 x 24mm Tubes
<b>Heat-On Multi-well System with Anodised Finish</b>	
RR62000	Heat-On Multi-well Block System Consisting of the following: 1 x RR62005 Heat-On Multi-well Holder 2 x RR62015 Heat-On 25ml Insert 2 x RR62020 Heat-On 50ml Insert 2 x RR62025 Heat-On 100ml Insert 1 x RR61085 Heat-On Adapter for 135mm Hotplates (IKA etc)
<b>Heat-On Multi-well Basic System with Anodised Finish</b>	
RR62001	Heat-On Multi-well Block System Consisting of the following: 1 x RR62005 Heat-On Multi-well Holder 1 x RR62015 Heat-On 25ml Insert 1 x RR62020 Heat-On 50ml Insert 1 x RR62025 Heat-On 100ml Insert 1 x RR61085 Heat-On Adapter for 135mm Hotplates (IKA etc)
<b>Heat-On Blocks with Anodised Finish</b>	
RR62031	Heat-On 100ml Block with Flask Sidearm Cutouts
RR63035	PTFE Safety Cover – for Heat-On 100ml Block
RR62032	Heat-On 200ml Florentine Block
RR62040	Heat-On 250ml Block
RR62041	Heat-On 250ml Block with Flask Sidearm Cutouts
RR62042	Heat-On 300ml Florentine Block
RR63040	PTFE Safety Cover - for Heat-On 200-300ml Block
RR62045	Heat-On 500ml Block
RR62046	Heat-On 500ml Florentine Block
RR63045	PTFE Safety Cover - for Heat-On 500ml Block
RR62050	Heat-On 1 Litre Block
RR62051	Heat-On 1 Litre Florentine Block
RR63050	PTFE Safety Cover - for Heat-On 1 Litre Block
RR62055	Heat-On 2 Litre Block
RR62056	Heat-On 2 Litre Florentine Block
RR62060	Heat-On 3 Litre Block
RR62061	Heat-On 4 Litre Block
RR62063	Heat-On 5 Litre Block

# Set-Up

## Select your stirring hotplate

Heat-On is compatible with stirring hotplates with a diameter of up to 145mm including the Radleys Carousel Stirring Hotplates and all popular Heidolph and IKA and stirring hotplates.

Carousel Standard

Carousel Tech

Carousel Advanced

## Do you need an adapter plate?

Heat-On blocks and Multi-well Holders have a single recess in the base to locate onto your stirrer hotplate. Some hotplates, such as the Carousel Stirring Hotplates and IKA RCT, have a smaller top plate diameter (135mm) so some movement may be noticed between the block and the top plate.

If you feel this is unacceptable an adapter plate can be purchased. This optional aluminium adapter plate for 135mm diameter top plates is available from Radleys:

Cat No	Description
RR61085	Heat-On Adapter Plate for 135mm Hotplates

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Always fully read the instructions of any third party equipment, such as stirring hotplates or temperature probes thoroughly, to ensure you are not contravening any safety recommendations or manufacturer's warranty when used in conjunction with Heat-On.

## Optimising the temperature performance of Heat-On

The heating performance of Heat-On will be affected by:

1. The limit of the maximum operating temperature of the stirring hotplate.
2. Heating power (W) of the stirring hotplate.
3. Flatness of the top plate on the stirring hotplate.
4. Ambient temperature.
5. Load to be heated.
6. Maximum housing temperature of the stirring hotplate.
7. Whether you use a separate temperature probe see overleaf.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Heat-On blocks can be used up to 200°C without damage to the block. Temperatures up to 260°C can be achieved but some minor surface degradation may occur.

### Warning

When placing oil baths, sand baths, polymer coated or aluminium blocks etc, onto a stirring hotplate the increased surface area will reflect heat back to the stirring hotplate. In such instances care should be taken to ensure that the hotplate's external housing does not exceed the maximum temperature recommended by the manufacturer. Such overheating can cause damage to internal and/or external parts and lead to failure of your stirring hotplate.

### Warning

Always refer to the manufacturer's operating instructions for your stirring hotplate before using it with Heat-On to be sure of any limitations or safety restrictions.

## Set-Up - continued

### Using a separate temperature probe (Pt1000)

For accurate temperature probe of your block and/or solution temperature we recommend using a temperature probe such as the RR91226 Pt1000 s/s Temperature Sensor.

All Heat-On blocks (plus Multi-well Holder) feature a hole (3.5mm ID) for inserting a temperature probe.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

#### Important Note

Please note that typically hotplates should be set between 5°C and 15°C above the boiling point of the solvent to achieve a gentle reflux.

#### Warning

Probe wire/cable should not come into contact with the Heat-On Block. Always rotate the block to ensure that the probe hole is nearest to the probe socket on the stirring hotplate.

### Optimising the stirring performance of your Heat-On

The Heat-On system uses the magnetic field of the stirring hotplate to stir your sample. Heat-On blocks have been optimised to work with leading brands of stirring hotplates with a top plate diameter of up to 145mm.

The stirring performance of Heat-On will be affected by:

1. The maximum stirring speed of the stirring hotplate.
2. Power and size of the magnet within the stirring hotplate.
3. Selection of an appropriate magnetic stirring bar for your chosen flask.
4. Viscosity of sample.

### Wide choice of magnetic stirring bars

The selection of an appropriate magnetic stirring bar for your chosen flask is key to the performance of stirring within the Heat-On block.

For more information on suitable stirring bars please contact Radleys or your local distributor.

RR98075	Cross Shape Stirring Bar 10mm, pk 40	RR95920	Oval Stirring Bar 40mm, pk 5
RR98091	Cross Shape Stirring Bar Rare Earth 16.5mm, pk 20	RR95921	Oval Stirring Bar 50mm, pk 5
RR95910	Cross Shape Stirring Bar 30mm, pk 5	RR95915	Hub Stirring Bar 45mm, pk 1
RR95905	Cylindrical Stirring Bar 8mm, pk 40	RR95916	Hub Stirring Bar 62mm, pk 1
RR98113	Cylindrical Stirring Bar Rare Earth 12mm, pk 20	RR95917	Hub Stirring Bar 70mm, pk 1
RR98070	Octagonal Stirring Bar 13mm, pk 20	RR95925	Spinring Stirring Bar 31mm, pk 1
RR99607	Pivot Ring Stirring Bar 12x4.5mm, pk 40	RR95926	Spinring Stirring Bar 44mm, pk 1
RR98071	Pivot Ring Stirring Bar 12x6mm, pk 40	RR95927	Spinring Stirring Bar 57mm, pk 1
RR99613	Pivot Ring Stirring Bar 15mm, pk 40	RR98095	PTFE Magnetic Stirring Bar Evaluation Kit, pk 30
RR98096	Elliptical Stirring Bar Rare Earth 10mm, pk 40	RR71200	PTFE Magnetic Stirring Bar Evaluation Kit - for flasks, pk 10
RR98097	Elliptical Stirring Bar Rare Earth 15mm, pk 20	RR98094	PTFE Magnetic Stirring Bar Retriever 350mm, pk 1
RR99064	Elliptical Stirring Bar Rare Earth 25mm, pk 10	RR98114	Magnetic Stirring Bar Restrainer, pk 1

## Set-Up - continued

### Beginning your experiment

In an ideal situation you should load your glassware into the block when the block is at ambient (room temperature).

However, in certain situations this may not be practical. In such circumstances you should note:

- You can use cold water to cool the block down rapidly, but take care with blocks above 100°C as steam is a potential hazard.
- It is generally safe to insert flasks into blocks that are below 50°C, but the flask should be removed before the block temperature drops below the start temperature as there is a small risk that the contraction of the block could crush the glass flask.
- In all eventualities flasks should be removed at or above the temperature they were originally inserted; e.g. in the case of a flask that is inserted at 80°C, it should be removed when or before the block cools to 80°C.


Heat-On blocks are available from 100ml to 5 litres and are designed to locate directly on to the stirring hotplate or adapter plate (see page 7).

- Select the appropriate Heat-On block for your chosen flask. Each Heat-On block is designed for a specific size of glassware and is identified accordingly with engraving on the side of the block.
- If using a 135mm hotplate, place adapter plate on hotplate.
- Place block onto hotplate/adapter plate, ensuring block is sitting level.
- If using a temperature probe rotate Heat-On block until probe hole is located at the rear of your stirring hotplate.

### Warning

Always refer to the manufacturer's operating instructions for your stirring hotplate before using it with Heat-On, to be sure of any limitations or safety restrictions.

### Important Note

Compatibility between Heat-On block and flask is essential for overall safety of this equipment (see page 12).

## Set-Up - continued

### Multi-well formats for smaller flasks and inserts

The unique Multi-well Holder is designed to hold either one or two inserts for flasks or tubes. The inserts are available for 10ml, 25ml, 50ml, 100ml or 150ml flasks as well as multi-tube inserts for 16mm  $\varnothing$ , 20mm  $\varnothing$  and 24mm  $\varnothing$  tubes or vials. Flask inserts feature cut-away sides to enable the use of flask with or without sidearms.

Simply place the Multi-well Holder directly on to the hotplate or adapter plate (see page 7) and select appropriate insert for your flask(s) or tube(s).

### Warning

Always refer to the manufacturer's operating instructions for your stirring hotplate before using it with Heat-On, to be sure of any limitations or safety restrictions.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Notes

- 1: Compatibility between Heat-On block and flask is essential for overall safety of this equipment (see page 12).
- 2: The Multi-well Holder only accepts ONE 150ml insert.

### Anodised blocks or fluoropolymer coated

The Heat-On blocks, holders and inserts can be purchased either polymer coated or anodised. The block's finish does not affect the temperature or stirring performance of the product, however the coatings do vary in regard to chemical resistance.

#### Anodised blocks

Anodised Heat-On blocks are resistant to the majority of solvents and splash resistant to dilute acids and alkalis at room temperature.

#### Fluoropolymer coated blocks

The fluoropolymer coated Heat-On blocks offer superior resistance to the majority of solvents, acids and alkalis.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Extended exposure to acids or alkalis will attack the surface of Heat-On. The stronger the concentration and the longer the exposure time the more chance and intensity of any attack. Heat will also speed up and intensify any attack; e.g., hot concentrated sulphuric acid will damage the fluoropolymer coated Heat-On, over a period of time.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

It is important to clean off any residual chemical spills immediately after they occur.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Heat-On blocks can be used up to 200°C without damage to the block. Temperatures up to 260°C can be achieved but some minor surface degradation may occur - the polymer coating may be damaged, and anodised blocks may become discoloured.

## Set-Up - continued

### Fitting optional safety lifting handles to Heat-On blocks and the Multi-well Holder

Safety lifting handles have a simple and secure quick release mechanism that allows the user to quickly attach the handles and move the block while it remains hot.

#### Warning

When lifting the block using the handles, always be sure to touch only the handles and not the block itself, which may be too hot to touch and may cause injury.

#### Warning

Always wear heat resistant gloves when handling hot blocks - even when using the lifting handles.

The handles are attached to the block or holder as follows:

1. Take handle and push button on top (as shown in picture to the right).
2. Slide handle into the stainless steel insert.
3. Release button.

To remove handle, simply press down on button and pull handle out.

#### Warning

Handles should always be removed from block when heating. If handles are left in block whilst heating they will become hot.

#### Warning

Extreme caution should always be used when lifting hot vessels or objects. Please check with your laboratory safety officer to ensure that you are complying with all relevant safety procedures.

**Many laboratories do not permit the handling of hot liquids or objects and therefore the use of the optional handles may not be permitted within such a laboratory.**

### Using a retort stand

To increase user safety Radleys recommend the use of a retort stand with Heat-On blocks, adding further stability to the system, particularly when using the larger blocks.

We find it best to locate the stand to the rear of your set-up thus avoiding any obstruction when using hotplate controls.

## Set-Up - continued

---

### Flask selection process

Heat-On is designed so that oversized flasks cannot fit through the slightly narrower opening at the top of the block well. This means that all suitable flasks should sit at the bottom of the block well, be a loose fit and as such spin freely in the Heat-On well.

When a flask is inserted into the well it is important to ensure that the flask touches the bottom of the well and does not 'catch or stick' on the sides of the block well. If the flask is 'catching or sticking' it may be deformed and may result in jamming and even cracking of the flask when the Heat-On block is cooling.

There is a simple rule when using glassware with the Heat-On: "If it does not spin, then throw it in the bin". If the flask does not spin freely then it is unsuitable for use in Heat-On and an alternative flask should be used.

During the design of Heat-On, Radleys evaluated and tested flasks from a wide range of glass manufacturers to ensure that Heat-On blocks will accept the widest possible range of sizes and styles. However, because the tolerances on mass produced glassware can vary greatly it is not possible to guarantee that all brands will fit. Indeed, experience shows us that glassware from a single manufacturer may vary considerably in size and shape from batch-to-batch and also from item-to-item. Customers may have to pre-select appropriate sized glassware to use with Heat-On.

---

### Warning

Factors that could lead to a flask failure:

#### Description of problem

1. The flask does not spin freely and therefore is deformed and is catching or sticking on the inside of the well. The flask may well crack when the Heat-On block cools.
2. There is grit or an object between the flask and the well.
3. The flask was sealed with a closure or cap and there was a build up of pressure in the flask during heating.
4. The flask was scratched, or cracked and failed due to weakness.
5. The flask was not properly annealed during manufacture and has cracked due to stress.

#### Remedy

1. *There is simple rule when using glassware with the Heat-on: "If it does not spin, then throw it in the bin". If the flask does not spin freely then it is unsuitable for use in Heat-On and an alternative flask should be used.*
2. *Always thoroughly clean the Heat-On block before and after use.*
3. *Never heat sealed glass vessels, unless they are specifically pressure rated.*
4. *Never use glassware that shows signs of scratching, cracking or has chips, as this causes inherent weaknesses in the glass strength.*
5. *Always buy glassware from reputable sources.*

## Set-Up - continued

### Set-up of PTFE Safety Covers

1. Follow the set-up procedure for the stirring hotplate and Heat-On block as described previously.
2. Ensure that you have the correct Safety Cover for the Heat-On block that you are using. You can check this by matching the size information inscribed on the PTFE Cover with that on the Heat-On block.
3. Carefully position the Heat-On Safety Cover over the Heat-On block. Check that the Cover is properly located and fitting on the block securely. There should be adequate clearance between the inner circumference of the Heat-On Safety Cover and the inner rim of the Heat-On bowl, such that there is no obstruction when glassware is positioned in the block. The PTFE cover should turn freely on the Heat-On block.
4. Rotate the cover until the access holes for the Heat-On safety lifting handles and temperature probe in the cover are correctly aligned with the respective holes in the Heat-On block.
5. Load glassware: Select the appropriate glassware for your application. Ensure that the glassware is fully compatible with the Heat-On block and Safety Cover to be used. Position the glassware in the Heat-On block. Check that the glassware is properly located in the block, and not constrained in anyway by the Heat-On Safety Cover.

### Set-up of PTFE Safety Covers for Multi-well Holder

A dedicated Multi-well Safety Cover is available for the Heat-On Multi-well Holder, which is compatible with the full range of inserts. Simply place the Multi-well holder directly on to the hotplate, position the PTFE Safety Cover and rotate so that the central cut out of the cover is correctly aligned with the locating rings in the Multi-well block holder and the access holes for Heat-On safety lifting handles and temperature probe.

## Stirring Hotplate Specifications and Selection Guide

Description	Radleys Carousel Stirring Hotplate Standard, Tech and Advanced	Heidolph Hei-Standard, Hei-Tec and Hei-End	IKA RCT Basic
Maximum Temp	300°C	300°C	310°C
Heating Power	800W	800W	600W
Max Stirring Speed	1400rpm	1400rpm	1200rpm
Top Plate Diameter	135mm	145mm	135mm
Is adapter plate recommended?	yes	no	yes

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Information is correct at time of publication. No responsibility can be taken for inaccuracy of stirring hotplate data. All tradenames are the property of their respective owners. No endorsement for the suitability of any stirring hotplate for use with the Heat-On blocks is given or implied by any of the respective stirring hotplate manufacturers.

# Warranty – Email Back

sales@radleys.co.uk

To qualify for your warranty please complete, scan and email this form to Radleys

Date of Purchase .....	
Supplier's Name and Address .....	
Product Batch/Serial No. (if shown) .....	
<b>Your Details</b>	
Mr   Mrs   Miss   Ms   Dr   Prof	
Name .....	
Position .....	
Dept .....	Building .....
Organisation .....	
Address 1 .....	
Address 2 .....	
Town/City .....	County/State .....
Country .....	Post/Zip Code .....
Telephone .....	Ext ..... Fax .....
Email .....	Website .....

**Type of Organisation; please tick all boxes relevant**

<input type="checkbox"/> Academic Institution	<input type="checkbox"/> Consumer Goods	<input type="checkbox"/> Defence/Military/Forensic	<input type="checkbox"/> Government	<input type="checkbox"/> Manufacturing/Industrial	<input type="checkbox"/> Polymers/Plastics
<input type="checkbox"/> Animal Health/Zoology	<input type="checkbox"/> Contract Lab	<input type="checkbox"/> Environmental/Water	<input type="checkbox"/> Hospital/Pharmacy	<input type="checkbox"/> Nuclear/Gas/Electric	<input type="checkbox"/> Process Engineering
<input type="checkbox"/> Agrochemical	<input type="checkbox"/> Contract Synthesis	<input type="checkbox"/> Flavours/Fragrances	<input type="checkbox"/> Instrum/Elect & Mech	<input type="checkbox"/> Petrochemical/Oil	<input type="checkbox"/> Research Institute
<input type="checkbox"/> Chemical Manufacture	<input type="checkbox"/> Cosmetics	<input type="checkbox"/> Food/Beverages	<input type="checkbox"/> Lab Equip Dealer/Mnf	<input type="checkbox"/> Pharma/Biotech/API	<input type="checkbox"/> Other.....

**Areas of Interest; please tick all boxes relevant**

<input type="checkbox"/> Analytical Chemistry	<input type="checkbox"/> Chromatography	<input type="checkbox"/> Estate & Facilities	<input type="checkbox"/> Health & Safety	<input type="checkbox"/> Organic Chemistry	<input type="checkbox"/> QC/QA
<input type="checkbox"/> Automation/HTS	<input type="checkbox"/> Clinical/Medical/Pathology	<input type="checkbox"/> Food & Agriculture	<input type="checkbox"/> Inorganic/Metallurgy	<input type="checkbox"/> Parallel Chem/Combi-Chem	<input type="checkbox"/> Sales & Marketing
<input type="checkbox"/> Biochemistry	<input type="checkbox"/> Construction	<input type="checkbox"/> Formulation	<input type="checkbox"/> Liquid Handling/MicroPlates	<input type="checkbox"/> Polymers & Oils	<input type="checkbox"/> Separation/SPE
<input type="checkbox"/> Biological Sciences	<input type="checkbox"/> Drug Discovery	<input type="checkbox"/> Geology	<input type="checkbox"/> Material Science	<input type="checkbox"/> Process Dev/Scale-up	<input type="checkbox"/> Support/Engineering
<input type="checkbox"/> Catalysis	<input type="checkbox"/> Environmental Health	<input type="checkbox"/> Glassblower	<input type="checkbox"/> Medical Devices	<input type="checkbox"/> Process Safety/Calorimetry	<input type="checkbox"/> Temperature Control
<input type="checkbox"/> Other.....			<input type="checkbox"/> Medicinal Chemistry	<input type="checkbox"/> Purchasing/Stores	<input type="checkbox"/> Veterinary

To request specific product information from Radleys please fill in below

<p><b>Benchtop and Hotplates</b></p> <input type="checkbox"/> Findenser Air Condenser <input type="checkbox"/> Heat-On Block System <input type="checkbox"/> Cool-It Insulated Bowls <input type="checkbox"/> StarFish Work Station <input type="checkbox"/> Carousel Stirring Hotplates <input type="checkbox"/> Overhead Stirrers <p><b>Jacketed Lab Reactors</b></p> <input type="checkbox"/> Reactor-Ready Lab Reactor <input type="checkbox"/> Reactor-Ready Duo Lab Reactor <input type="checkbox"/> Reactor-Ready Pilot Lab Reactor <input type="checkbox"/> Custom Jacketed Reaction Systems	<p><b>Parallel Reaction Stations</b></p> <input type="checkbox"/> Carousel 12 Plus Reaction Station <input type="checkbox"/> Cooled Carousel 12 Reaction Station <input type="checkbox"/> Carousel 6 Plus Reaction System <input type="checkbox"/> Cooled Carousel 6 Plus Reaction Station <input type="checkbox"/> Carousel Work-Up Station <input type="checkbox"/> GreenHouse Plus Parallel Synthesiser <input type="checkbox"/> GreenHouse Work-Up Station <input type="checkbox"/> GreenHouse Blowdown Evaporator <input type="checkbox"/> Tomado Overhead Stirring System <input type="checkbox"/> Breeze Heating/Cooling Work Station <input type="checkbox"/> Storm Heating/Cooling Work Station	<p><b>Software</b></p> <input type="checkbox"/> AVA Lab Control Software <input type="checkbox"/> Level 1/2 <input type="checkbox"/> Level 3/4 <input type="checkbox"/> Data Hub <p><b>Automated Reaction Stations</b></p> <input type="checkbox"/> Mya 4 Reaction Station <p><b>Other</b></p> <input type="checkbox"/> Huber..... <input type="checkbox"/> Heidolph..... <input type="checkbox"/> Other.....
--	--	--

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**RB Radleys LTD**