

VIBRATORY DISC MILL RS 200



The Vibratory Disc Mill RS 200 is suitable for the **extremely quick, loss-free and reproducible grinding** of medium-hard, brittle and fibrous materials to analytical fineness. The instrument runs **steadily and smoothly**, even with heavy grinding sets, at maximum speed thanks to the new Stabilized-Plane-Drive.

With its **robust design**, the RS 200 has proven to be ideal for use in the building materials sector (cement), in geology, mineralogy, metallurgy and in power plants. Due to the **high end fineness** and speed the RS 200 is the perfect mill when it comes to **preparing samples for spectral analysis**.

APPLICATION EXAMPLES

cement, cement clinker, ceramics, coal, coke, concrete, corundum, electronic components, glass, metal oxides, minerals, ores, plant materials, silicates, slag, soils, ...

PRODUCT ADVANTAGES

- | reproducible results due to "Stabilized-Plane-Drive" (avoids tumbling of the grinding set)
- | variable speed 700 - 1,500 min⁻¹
- | extremely short grinding time
- | agate and tungsten carbide recognition for automatic speed reduction
- | 10 SOPs can be stored
- | no sample loss thanks to optimum o-ring seal
- | wide range of materials for contamination free grinding
- | noise-insulated grinding chamber
- | new quick-action grinding set clamping device for convenient and safe tightening
- | safety interlock
- | easy 1-button operation with colored graphics display
- | new carry handle
- | new ergonomic design, heavy grinding sets can be inserted closer to the user's body and slid into the machine, leading to less strenuous laboratory work especially for high throughput users

FEATURES


Applications	size reduction, mixing, tritulating
Field of application	construction materials, environment / recycling, geology / metallurgy, glass / ceramics
Feed material	medium-hard, hard, brittle, fibrous
Size reduction principle	pressure, friction
Material feed size*	< 15 mm
Final fineness*	< 20 µm
Batch size / feed quantity*	15- 250 ml depending on size of grinding set
Speed at 50 Hz (60 Hz)	700 - 1,500 min ⁻¹ , continuously adjustable
Material of grinding tools	hardened steel, tungsten carbide, agate, zirconium oxide, steel 1.1740 (for heavy-metal free grinding)
Grinding jar sizes	50 ml / 100 ml / 250 ml
Setting of grinding time	digital, 00:01 to 99:59
Storable SOPs	10
Drive	3-phase asynchronous motor with frequency converter
Drive power	1.5 kW
Electrical supply data	different voltages
Power connection	1-phase
Protection code	IP 40
W x H x D closed	836 x 1220 x 780 mm
Net weight	~ 210 kg (without grinding set)
Standards	CE

FUNCTION PRINCIPLE

The Vibratory Disc Mill RS 200 comminutes by **impact and friction**. The grinding set is firmly attached to the vibration plate with a quick-action lever. The plate with the grinding set is subjected to circular horizontal vibrations.

The centrifugal force acting on the grinding rings in the dish results in extreme pressure, impact and frictional forces acting on the sample, producing analytical fineness in 1-3 minutes. The circular vibrations are produced by a frequency controlled 1.5 kW 3-phase motor.

A sensor recognizes the presence of agate grinding sets and automatically limits the speed to 700 min⁻¹ to avoid damaging the agate. The cover of the noiseinsulated grinding chamber has a safety interlock and can only be opened when the mill is at a standstill.



[Click to view video](#)



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.