

STEAM INDICATOR (Autoclave) TAPE For Monitoring Steam Sterilisation Processes

Rodwell Part No. **AC985 (885-0306)**
 OEM Code: TAP-SCS-E



Product Description

Excelsior Steam autoclave indicator tape contains no lead or other toxic heavy metals. The autoclave tape is manufactured for closing packages, instrument wraps, and other items requiring steam sterilisation. The autoclave tape can be utilised for monitoring steam sterilisation processes (pre-vacuum and gravity) up to 140°C. The tape adhesive can adhere to paper, plastic, glass, metal, and non-woven materials.

Physical Properties

Process	Steam
Tape Dimensions	19mm x 55m (0.75" x 60 yards) Roll Outer Diameter 114mm (4.5")
Packaging	48 Rolls
Chemical Indicator	Initial Colour: Green Signal Colour: Orange Chemical Indicator Inks contain no lead or other toxic heavy metals

Indications for Use

Type 4 Multi-critical Process Variable Indicator :

- 121°C, ≥ 15 minutes (pre-vacuum or gravity)
- 132°C, 3.5 minutes (pre-vacuum)
- 134°C, 3 minutes (pre-vacuum)

Instructions for Use

Use a minimum of a three (3) inch section of tape to close packages/items intended for steam sterilisation. Process packages/items as instructed in the autoclave/steriliser validation or manual.

Upon exposure to high temperature steam, the chemical indicator will transition from green to orange. The transition colour may vary depending on the load configuration, length and conditions of exposure. A colour transition from green to a shade of orange or tan provides indication of exposure to high temperature steam.

The chemical reaction which causes the colour transition is a steam specific reaction and is irreversible.



*Colours shown are representations of printed ink initial and signal colours but may vary from actual use.









The signal colour achieved from exposure to steam may vary from the example above due to differences in processing parameters (i.e. load content, steam quality, cycle time, temperature, etc.). For a Type 1 Process Indicator, any colour change produced during exposure to steam which is different from the initial colour is considered acceptable.

Compliance

ISO 11140-1:2014 Sterilization of health care products – Chemical Indicators- Part 1: General Requirements for Type 1 Chemical Process Indicators & Type 4 Multi-critical Process Variable Indicators.

Storage and Shelf Life

	15°C to 30°C		Keep away from sunlight
	10% to 80% relative humidity		Keep Dry
Shelf Life	4 years from the date of manufacture. The date of manufacture is based on the day the indicating ink is applied to the substrate. The remaining shelf life upon receipt will be shorter than 4 years		Do not reuse
	Do not use the product if the indicator has transitioned from green to orange. Keep away from sterilants. Do not use after expiration date.		

Disposal

Discard as general waste.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.