



Battery: LR-3

PRODUCT DESCRIPTION

A combined temperature and humidity meter ideal for environmental monitoring of clean rooms, laboratories, computer rooms and storage rooms. The large LCD display is easy to read at a distance.

- -10 to 50°C x 0.1°
- 30 to 90%RH x 1%
- Displays temperature in °C or °F
- Simultaneous display of humidity and temperature
- 5-Level Comfort Indicator
- 75 x 63mm large digit LCD display
- Min/max memory for both temperature and humidity
- Reset clears memory for both temperature and humidity
- Impact-resistant ABS case with a flip out desk stand for free standing applications and key-way for wall mounting
- Supplied with battery

TECHNICAL SPECIFICATIONS

	Temperature	Humidity
Range	-10 to +50°C	30 to 90%RH
Resolution	0.1°C	1% RH
Accuracy	±1°C (0 to 40°C) ±2°C (others)	±5%RH (50-80%RH) ±8%RH (others)
Power	1 x 1.5V AAA alkaline battery (LR-3)	
Dimensions	100 x 20 x 135mm	
Weight	143g	



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.