



# SL-CLASSIC-UB ◀ 140 cm ▶



**Safe and approved storage of aggressive, non flammable hazardous materials in working areas**

### Function / construction:

- **robust construction and longevity:** made out of highly resistant melamin resin special plates
- **user-friendly:** access to all containers, 135° opening of the wing door, interior of the cabinet completely visible
- **no unauthorised use:** doors lockable with cylinder lock and locking state indicator (red/green)
- **mobile:** castors with plinth (optionally available), easy moving of the cabinet
- **ventilation:** integrated air ducts ready for connection (DN 50) to a technical exhaust system

### Available equipment

- pull-out shelves with removable plastic sump
- **model SL.060.140.UB:** right compartment with integral shelf for material storage

#### Cabinet for acids and alkalis

##### SL-CLASSIC-UB model SL.060.140.UB.3T

in Light grey (RAL 7035), interior equipment with 4 x pull-out shelf (panel/PP)

Order No. 30616-001-30618



#### Cabinet for acids and alkalis

##### SL-CLASSIC-UB model SL.060.059.UB.T

in Light grey (RAL 7035), interior equipment with 2 x pull-out shelf (panel/PP)

Order No. 30622-001-30625



#### Cabinet for acids and alkalis

##### SL-CLASSIC-UB model SL.060.110.UB.2T


in Light grey (RAL 7035), interior equipment with 4 x pull-out shelf (panel/PP)

Order No. 30619-001-30621


Technical data		SL.060.140.UB.3T	SL.060.110.UB.2T	SL.060.059.UB.T
External dimensions W x D x H	mm	1400 x 570 x 600	1100 x 570 x 600	590 x 570 x 600
Internal dimensions W x D x H	mm	left: 489 x 531 x 562 central: 489 x 531 x 562 right: 281 x 531 x 562	all compartments: 489 x 531 x 562	488 x 531 x 562
Weight without interior equip.	kg	80	60	40
Extraction air	DN	50	50	50
Air change <b>30 times</b>	m <sup>3</sup> /h	9	9	4
Diff. pressure (air change <b>30 times</b> )	Pa	10	10	5
Distributed load	kg/m <sup>2</sup>	82	93	151
Maximum load	kg	100	100	100




**SL.060.140.UB.3T Front view**




**Side view**




**Top view**




**SL.060.110.UB.2T Front view**




**Side view**




**Top view**




**SL.060.059.UB.T Front view**



**Side view**



**Top view**



SL-LINE

Model	Information / equipment	Order No.
SL.060.140.UB.3T	without interior equipment, Underbench cabinet, with 3 wing doors	30616
SL.060.110.UB.2T	without interior equipment, Underbench cabinet, with 2 wing doors	30619
SL.060.059.UB.T	without interior equipment, Underbench cabinet, left hinged, with wing door	30622
SL.060.059.UB.TR	without interior equipment, Underbench cabinet, right hinged, with wing door	30623

Body colour	Door colour	Order No.
Light grey RAL 7035	Light grey RAL 7035	001 <b>4</b>

Interior equipment packages (can only be ordered in combination with the cabinet) Width 1400 mm	Material	Order No.
left: 1x pull-out shelf (V=11.0L) right: 1x pull-out shelf (V=11.0L)	panel/PP	30617 <b>2</b>
left: 2x pull-out shelf (V=11.0L) right: 2x pull-out shelf (V=11.0L)	panel/PP	30618 <b>2</b>

Interior equipment packages (can only be ordered in combination with the cabinet) Width 1100 mm	Material	Order No.
left: 1x pull-out shelf (V=11.0L) right: 1x pull-out shelf (V=11.0L)	panel/PP	30620 <b>2</b>
left: 2x pull-out shelf (V=11.0L) right: 2x pull-out shelf (V=11.0L)	panel/PP	30621 <b>2</b>

Interior equipment packages (can only be ordered in combination with the cabinet) Width 590 mm	Material	Order No.
1x pull-out shelf (V=11.0L)	panel/PP	30624 <b>2</b>
2x pull-out shelf (V=11.0L)	panel/PP	30625 <b>2</b>

Interior equipment	Material	For cabinet width	1)	2)	3)	Order No.
pull-out shelf, with PP sump	panel/PP grey	590, 1100 and 1400	25.00	18.00	16.36	30770 <b>2</b>
castors with plinth, height = 30 mm	sheet steel powder-coated RAL 7035	1400				30374 <b>2</b>
castors with plinth, height = 30 mm	sheet steel powder-coated RAL 7035	1100				30362 <b>2</b>
castors with plinth, height = 30 mm	sheet steel powder-coated RAL 7035	590				30478 <b>2</b>

Accessories	Order No.
extraction module without exhaust air monitoring for wall mounting	17178 <b>2</b>

\* usable storage space W x D (mm)  
1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.