



# DYNAREAD Microplate Absorbance Reader



## Key benefits of the **DYNAREAD**

- Absorbance reading and data analysis of 96-well microplates for ELISA
- Linear shaking function, with adjustable speed and amplitude
- Reader control and results data analysis with connected computer system
- Compatible with well-established and widely used Dynex Revelation software
- Fully comprehensive service and support available
- Accommodates up to 8 optical filters (4 filters included)
- Wide dynamic range of 0.000 – 4.000 OD
- Long-life, reliable and stable LED light source
- End-point and kinetic read methods
- High-speed 12-channel reading
- Wavelength range 405 – 790 nm
- CE IVD marked





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.