

MORTEXER™

MORE Power, MORE Flexibility... MORE Vortexer!



- The only vortexer available with **Multi-Head™**
- **Powerful Motor for INSTANT Vortexing**
- **Q-Drive™ Dynamic Balance System**
- **Touch or continuous operation**



Unique **Multi-Head** allows for standard vortexing **PLUS** holds up to 8 microtubes

The *Mortexer™* is MORE than your standard vortex mixer. The unique design of the Multi-Head™ (included) features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes.

The drive system is balanced with Q-Drive™ technology, providing a smooth vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation. Increased speed, up to 3400 rpm produces an instant vortex of even the largest sample sizes, including full 50ml tubes.

Other features include: "no-walk" suction feet, cold room/incubator compatibility and spill-proof electronic controls. Two year warranty.

Technical Data:

Speed Range:	200 to 3200 rpm
Operating Modes:	Touch or continuous
Operating Temp. Range:	+4°C to +65°C
Dimensions:	13 x 16 x 17 cm
	5 x 6.3 x 6.7 in.
Weight:	3.8 kg / 8.4 lb
Electrical:	115V AC, 60 Hz, 1.5A
	230V AC, 50 Hz, 0.75A
	230V AC, 60 Hz, 0.75A
Warranty:	2 years



Ordering Information:

BV1005* *Mortexer™* Vortex Mixer (with Multi-Head)

Optional Accessories:

BV1000-FLAT	Optional 3 inch flat head
BV1000-COMBO	Optional head for tubes and microplate. 66 tubes (38 x 1.5ml and 28 x 0.5ml) or 1 microplate
BV1000-H15	Optional head for 12 x 1.5/2.0ml, horizontal
BV1000-H150	Optional head for 4 x 15ml, horizontal
BV1000-H500	Optional head for 2 x 50ml, horizontal

* To order a product in 230V / 50Hz, please add -E to item number.
To order a product in 230V / 60Hz, please add -E2 to item number.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.