


EXPellPLUS

LOW RETENTION TIPS

Thanks to the unique composition of the plastic polymer, **CAPP's** premium **EXPellPlus** low retention pipette tips have an ultra-smooth inner surface, minimising any sample binding. The polymer composition, along with the high quality Swiss molds on which the tips are manufactured, enables **CAPP** to obtain the low retention properties without the use of any surface additives.

EXPellPlus low retention tips are ideal for pipetting viscous samples, which are normally difficult to dispense in full. Samples such as plasma, enzymes, foaming liquids, detergents, TRIS or BSA can be handled comfortably and without any waste with **CAPP** low retention tips. Additionally, **EXPellPlus** low retention tips can significantly reduce the waste of expensive reagents. **EXPellPlus** low retention tips are the perfect choice for volumes below 1mL and help to improve the accuracy and reproducibility of pipetting results. Even the smallest amount of liquid remaining in the tip after dispensing can have a negative impact on the pipetting results. **CAPP EXPellPlus** low retention tips are manufactured in a clean room facility and are all RNase, DNase and Pyrogen free certified and sterile.

Below picture shows a comparison of the liquid remaining after dispensing in a standard tip and in an **EXPellPlus** low retention tip. **EXPellPlus** low retention tips are universally compatible with the majority of other pipette brands.



Ordering Information

EXPellPlus Low Retention Non-Filter Tips

Cat. No.	Description
5030010C	ExpellPlus 10µL, pre-sterile, clear, bag, 20x1000 pcs.
5030012C	ExpellPlus 10µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped
5030020C	ExpellPlus 10µL, pre-sterile, hinged racks, 5x10x96 pcs.
5030040C	ExpellPlus 10µL XL, pre-sterile, clear, bag, 20x1000 pcs.
5030042C	ExpellPlus 10µL XL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped
5030050C	ExpellPlus 10µL XL, pre-sterile, hinged racks, 5x10x96 pcs.
5030051C	ExpellPlus 10ul XL, extra narrow, pre-sterile, hinged racks, 5x10x96 pcs.
5030070C	ExpellPlus 200µL, pre-sterile, clear, bag, 30x1000 pcs.
5030072C	ExpellPlus 200µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped
5030080C	ExpellPlus 200µL, pre-sterile, hinged racks, 5x10x96 pcs.
5130110C	ExpellPlus 300µL, pre-sterile, hinged racks, 5x10x96 pcs.
5030130C	ExpellPlus 1000µl (1250µL), pre-sterile, clear, bag, 10x1000 pcs.

EXPellPlus Low Retention Filter Tips

Cat. No.	Description
5030030C	ExpellPlus 10µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030060C	ExpellPlus 10µL XL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030061C	ExpellPlus 10µL XL, extra narrow, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030062C	ExpellPlus 20µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030066C	ExpellPlus 100µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030090C	ExpellPlus 200µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030115C	ExpellPlus 300µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs.
5030150C	ExpellPlus 1000µL, pre-sterile w/ filter, hinged racks, 4x8x96 pcs.



COMPATIBILITY CHART

Description	Cat. No.	CAPP			Gilson			Eppendorf			Brand			Thermo			Sartorius			HTL		
		Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic	Pipetman / Pipetman G	Research Plus	Transfipette, Transfipette S	Fimmpipette E1 / Novus	Proline Plus / Picus	Labmate Pro (same as Labnet BioPette)	Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic	Pipetman / Pipetman G	Research Plus	Transfipette, Transfipette S	Fimmpipette E1 / Novus	Proline Plus / Picus	Labmate Pro (same as Labnet BioPette)	Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic	Pipetman / Pipetman G	Research Plus	Transfipette, Transfipette S	Fimmpipette E1 / Novus	Proline Plus / Picus	Labmate Pro (same as Labnet BioPette)
CAPP Expell Filter Tips																						
Expell 10 µL, w/ filter, clear, sterile	5130030C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10 µL, w/ filter, clear	5131030C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10 µL XL, w/ filter, clear, sterile	5130060C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10 µL XL, w/ filter, clear	5131060C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10/20 µL XL, extra narrow, w/ filter, clear, sterile	5130061C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10/20 µL XL, extra narrow, w/ filter, clear	5131061C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 20 µL, w/ filter, clear, sterile	5130062C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 20 µL, w/ filter, clear	5131062C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 100 µL, w/ filter, clear, sterile	5130066C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 100 µL, w/ filter, clear	5131066C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 200 µL, w/ filter, clear, sterile	5130090C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 200 µL, w/ filter, clear	5131090C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1000 µL, w/ filter, clear, sterile	5130150C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1000 µL, w/ filter, clear	5131150C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1250 µL XL, w/ filter, clear, sterile	5130123C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
CAPP Expell and ExpellPlus individually wrapped Pipette Tips																						
ExpellPlus 10 µL, clear, sterile, indiv. wrapped	5030012C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 10 µL XL, clear, sterile, indiv. wrapped	5030042C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 200 µL, clear, sterile, indiv. wrapped	5030072C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1000 µL (1250 µL), clear, sterile, indiv. wrapped	5130132C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
CAPP Expell and ExpellPlus Filter Tips in PaperBox Packs																						
ExpellPlus 10 µL, w/ filter, clear, sterile	5030029C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile	5030059C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 20 µL, w/ filter, clear, sterile	5030063C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 100 µL, w/ filter, clear, sterile	5030064C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 200 µL, w/ filter, clear, sterile	5030089C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 300 µL, w/ filter, clear, sterile	5030114C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1000 µL, w/ filter, clear, sterile	5130149C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ExpellPlus 1250 µL XL, w/ filter, clear, sterile	5130119C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
CAPP Expell Pipette Tips Reload System																						
Expell 10 µL, clear, reload	5030035C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 10 µL XL, extra narrow, clear, reload	5030055C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 200 µL, clear, reload	5030095C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 300 µL, clear, reload	5130111C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Expell 1000 µL (1250 µL), clear, reload	5130155C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■

This compatibility chart is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes other than those of CAPP. Specifications like inaccuracy and imprecision values have not been checked. For a complete compatibility list please contact info@cappahn.com.

■ compatible

■ X not compatible





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.