

# CAPP Solo

## SINGLE CHANNEL PIPETTES

**CAPP Solo** single channel pipettes are cost effective and robust liquid-handling tools designed for academic, diagnostic and start-up laboratories with limited budgets.

These pipettes meet industry standards of accuracy and precision, ensuring no compromises in all liquid handling protocols. The single channel pipettes are universally compatible with **CAPP EXPellPlus** pipette tips and most common tip brands on the market.

### Features:

- Colour coded for fast and easy visual identification of pipettes with different volume capacities
- Fine-tuned to ensure high levels of pipetting accuracy, precision, and long-term calibration retention
- Exceptionally light pipette chassis reduces the risk of repetitive strain injuries
  - Lockable counter-mechanism prevents accidental turning while pipetting
  - Toughened to withstand the toughest laboratory conditions including chemical exposure, drops, and autoclaving amongst others
  - Tip cone designed for universal tip compatibility
  - Backed by a two-year warranty



## Ordering Information

### CAPP Solo Fixed Volume

Cat. No.	Description
S5-1F	CAPPSolo pipette, fixed vol. 5µL
S10-1F	CAPPSolo pipette, fixed vol. 10µL
S20-1F	CAPPSolo pipette, fixed vol. 20µL
S25-1F	CAPPSolo pipette, fixed vol. 25µL
S50-1F	CAPPSolo pipette, fixed vol. 50µL
S100-1F	CAPPSolo pipette, fixed vol. 100µL
S200-1F	CAPPSolo pipette, fixed vol. 200µL
S250-1F	CAPPSolo pipette, fixed vol. 250µL
S500-1F	CAPPSolo pipette, fixed vol. 500µL
S1000-1F	CAPPSolo pipette, fixed vol. 1000µL
S2000-1F	CAPPSolo pipette, fixed vol. 2 mL
S2500-1F	CAPPSolo pipette, fixed vol. 2,5 mL
S5000-1F	CAPPSolo pipette, fixed vol. 5 mL
S10000-1F	CAPPSolo pipette, fixed vol. 10 mL

### CAPP Solo Variable Volume

Cat. No.	Description
S2-1	CAPPSolo pipette, variable vol. 0.1 – 2.5µL
S10-1	CAPPSolo pipette, variable vol. 0.5 – 10µL
S20-1	CAPPSolo pipette, variable vol. 2 – 20µL
S50-1	CAPPSolo pipette, variable vol. 5 – 50µL
S100-1	CAPPSolo pipette, variable vol. 10 – 100µL
S200-1	CAPPSolo pipette, variable vol. 20 – 200µL
S1000-1	CAPPSolo pipette, variable vol. 100 – 1000µL
S5000-1	CAPPSolo pipette, variable vol. 0.5 – 5 mL
S10000-1	CAPPSolo pipette variable vol. 1 – 10 mL





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.