

CAPP Rondo

SINGLE PLATE SHAKER



The **CAPP Rondo Plate Shaker** is an ideal choice of microplate shaker for any test environment. With a tiny footprint, as small as the 96-plate size, the shaker is microprocessor-controlled for higher accuracy. The in-built counterbalance ensures a stable, quiet and long life operation.

Highlights:

- Comes with interchangeable clamps that provide universal fitting for different microplates
- Provided with an accessory tube platform to accommodate 20 microtubes of 1.5/2mL
- With a wide speed setting range from 300 to 1.800 RPM and orbital motion of 2mm (1mm radius), the **CAPP Rondo Plate Shaker** can perform numerous applications including PCR, ELISA, enzyme immunoassays, protein synthesis and pharmaceutical profiling
- Vortexing motion thoroughly mixes each well
- Can be used inside cold rooms or fume hoods (up to 40°C)

Specifications:

| | |
|-------------------------|---|
| Speed Setting | 300 to 1.800 RPM |
| Digital Display | Yes |
| Shaking radius | 1mm |
| Motor Type | Brushless DC |
| Load Capacity | 2kg |
| In-built Counterbalance | Yes |
| Power Supply | Universal (100-240V, 50/60 Hz.) |
| Dimensions (W x D x H) | 127 x 85 x 86mm |
| Weight | 1.36kg |
| Accessories | Replaceable clamps to use different types of standard microplates |

Ordering Information

| Cat. No. | Description |
|----------|--|
| CRP-18X | CAPPRondo plate shaker, 300 to 1.800 RPM |





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.