


# identitapes

Self-adhesive IdentiTapes are designed to provide the perfect solution to a wide variety of labelling requirements


Cleaver IdentiTapes are oil-proof, waterproof and acid resistant tapes which will withstand temperature extremes of -20°C to over +121°C. They will adhere to surfaces such as glass, plastic, metal, paper, rubber and any other dry surface yet can be easily removed with no adhesive residue being left on the original surface. IdentiTapes can be written on with ball point pen, pencil, felt markers.

## KEY FEATURES

- Identify samples easily
- Oil, water and acid resistant
- Paper based
- Withstands extremes in temperatures, -20°C to +121.1°C
- Ideal for use in a freezer, refrigerator or incubator
- No adhesive residue left when removed
- Surface can be written on with pen, pencil or lab marker




Sharp cutter bar



Acrylic Tape Dispenser, TD1

### 55 m Core Adapters, CA55

IdentiTape Dispensers are designed to accommodate tapes of both core diameters. The central spindle takes the smaller 12.7 m rolls directly while by using the supplied adapters, the larger 55 m rolls, shown in the adjacent photograph, allows both



# coolcube

CoolCube is designed to store samples safely on the lab bench without sample degradation caused by temperature increase or fluctuation.

Simply place it in a freezer overnight and it will be ready to use.

CoolCube can be used as a convenient PCR workstation for 12 standard 1.5/2.0ml microtubes with one 96 well PCR plate, or with 12x 0.2ml strips or 96x0.2ml tubes. The other side may be used to hold up to 36 standard 1.5/2.0ml microtubes.



## KEY FEATURES

- Maintains 0°C for 4 hours
- Holds both microtubes and a PCR plate
- Can be used as a PCR workstation

# ice baths

IceCool Baths are constructed of durable plastic and offer both robustness and long life. They will also withstand reasonable use with dry ice. Lids are available for the large bucket.



## KEY FEATURES

- Durable plastic
- Easily cleaned
- Excellent long term thermal storage

## ORDERING INFORMATION

<b>identitapes</b>									
<b>Length</b>	<b>Width</b>	WHITE	GREEN	YELLOW	BLUE	RED	ORANGE	PINK	VIOLET
12.7m	12.7mm	ST-12/W	ST-12/G	ST-12/Y	ST-12/B	ST-12/R	ST-12/O	ST-12/P	ST-12/V
12.7m	19.0mm	ST-34/W	ST-34/G	ST-34/Y	ST-34/B	ST-34/R	ST-34/O	ST-34/P	ST-34/V
12.7m	25.4mm	ST-10/W	ST-10/G	ST-10/Y	ST-10/B	ST-10/R	ST-10/O	ST-10/P	ST-10/V
55m	12.7mm	ST-120/W	ST-120/G	ST-120/Y	ST-120/B	ST-120/R	ST-120/O	ST-120/P	ST-120/V
55m	19.0mm	ST-340/W	ST-340/G	ST-340/Y	ST-340/B	ST-340/R	ST-340/O	ST-340/P	ST-340/V
55m	25.4mm	ST-100/W	ST-100/G	ST-100/Y	ST-100/B	ST-100/R	ST-100/O	ST-100/P	ST-100/V
<b>TD1</b>	<b>Acrylic Tape Dispenser, including CA55 core adapters</b>					<b>CA55</b>	<b>Additional core adapters, pack/5</b>		
<b>CSLIB</b>	Large Ice Bucket, 4L - White					<b>CSMIB</b>	Mini Ice Bucket, 1.4L - White		
<b>CSLIBLID</b>	Lid For Large Ice Bucket - White					<b>CSL-COOLCUBE</b>	Coolcube Storage Rack		



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.