



## MINIVORTEX

---

### eCommerce

#### Ecommerce Description

The Mini Vortexer's small footprint with maximum power offers immediate vortex shaking for a variety of tube sizes. A simple, touch start operation enables exceptional mixing performance using one hand.

- ABS and silicone materials to resist most commonly used chemicals
- Low profile design and small footprint (only 4.3" wide) save valuable space on a bench top and in fume hoods
- Rapid vortex action saves time and enables quick mixing of samples
- Modern design allows easy one handed touch start operation and shaking stops when the tube is lifted
- When performing repetitive procedures, a low activation force helps prevent user fatigue
- Maintenance free design and 3 year warranty helps reduce the cost of ownership
- Stable platform reduces the chance of dropping a sample during operation
- Quiet operation and low vibration helps users stay focused

Package Includes: 1 Mini Vortexer, 2 Silicone Caps, 1 Low voltage, double insulated power adapter with 4 interchangeable plugs

---

#### Specifications Table

Specifications: Speed: 4500 rpm

Orbit: 6 mm

Capacity: Up to 50 mL tube

Noise level: Full speed, unloaded: 53 dB

Unit Power: Input: 12 VDC, approx. 7.8 W (under load)

Adapter:

Input: 100-240 VAC, 50/60 Hz, 0.3 A

Output: 12 VDC, 1A with 4 interchangeable plug adapters. For indoor use, DoubleInsulated, Level VI efficiency,

Recognized marks: UL (Canada and U.S.), GS, FCC, CE, and C-Tick

Certifications: CE, CB Scheme, and SGS (Canada and U.S.)

---

### Ecommerce Short Description

- ABS and silicone materials to resist most commonly used chemicals
- Low profile design and small footprint (only 4.3" wide) save valuable space on a bench top and in fume hoods
- Rapid vortex action saves time and enables quick mixing of samples
- Modern design allows easy one handed touch start operation and shaking stops when the tube is lifted
- When performing repetitive procedures, a low activation force helps prevent user fatigue
- Maintenance free design and 3 year warranty helps reduce the cost of ownership
- Stable platform reduces the chance of dropping a sample during operation
- Quiet operation and low vibration helps users stay focused

---

### Product Includes

1 Mini Vortexer, 2 Silicone Caps, 1 Low voltage, double insulated power adapter with 4 interchangeable plugs

---

### URL

<https://www.cleaverscientific.com/electrophoresis-products/mini-vortex-mixer/>

---

## Product Management

### GTIN

5055323259982

---

### Product Group

Shakers and Vortexes

---

## General

### Unit Dimensions (cm)

12.954 x 10.922 x 7.62

---

### Capacity

2mL to 50 mL tube

---

### Weight (Kg)

0.90

---

### Construction

ABS & Silicone

---

## Files

Web Optimised Images

---

## Shipping Information

Shipping Length (cm)	20.32
----------------------	-------

---

UNSPC Code	41103814
------------	----------

---

Shipping Height (cm)	8.89
----------------------	------

---

Shipping Weight (Kg)	1.25
----------------------	------

---

Shipping Width (cm)	15.24
---------------------	-------

---

## Electrical

Power	7.8W
-------	------

---

Max current consumption	1A
-------------------------	----

---

Input Voltage	12VDC/ 110-240 VAC
---------------	--------------------

---

Compliance Marks	CE, SGS, RoHS2
------------------	----------------



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

