

## MXS-10 Roller Mixer

*Clifton Range*<sup>®</sup>



Heavy duty Clifton Roller Mixer for vials, small bottles, and test tubes. Durable design for mixing solid and liquid suspensions and viscous liquids, featuring rugged sprocket and chain drive, coupled with ball race roller mountings for smooth trouble free operation and can run continuously for many hours unattended using powerful DC motor.

Ten coated aluminium rollers run parallel to create a smooth rolling action between 10 and 60 rpm.

### MXS-10

GENERAL	
Dimensions	400w x 500d x 120h mm
Mixing Type	Rolling, no tilt.
Mixing Capacity	10 Rollers, Each pair of rollers can mix up to 1kg in weight
Mixing Speed	10 - 60 rpm, analogue control
Roller Size	30mm diameter x 300mm long
Permissible Mixing Time	100% usage
Ambient Temp Range	0 - 60°C
Stacking	Roller Mixers can be stacked 2 units high (accessory kit required)
Weight	7kg
Voltage	100 - 230v 50/60Hz





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.