

## 1000 Series Boiling Baths



1026ER



1026E



Round stainless steel electrically heated, the 1026ER regulates from a gentle simmer to a vigorous boil, whereas the 1026E has a timed temperature. Fitted with concentric ringed lid to vary the hole size and fitted with a sheathed element and a resettable thermal cut-out for safety.

Features an integral constant level device to maintain constant water level during boiling and supplied with a set of rings to vary the diameter of aperture to accommodate different sizes of vessels e.g. evaporating dishes and glassware, including round-bottom flasks.

### 1026E

### 1026ER

| GENERAL                 | 1026E                       | 1026ER  |
|-------------------------|-----------------------------|---|
| Temp. Range             | Boiling at 100°C            | Energy regulator controls the amount of boiling from a gentle simmer to vigorous boil |
| Temp. Control           | Constant boiling no control | Simmer to boiling at 100°C  |
| Heater Power (220-240V) | 1000W                       |   |
| Voltage                 | 230V models 50-60Hz         |   |
| DIMENSIONS              |                             |   |
| Tank Capacity           | 3.24 litres                 |   |
| External (mm)           | 200w x 130d                 |   |

By removing concentric rings, the hole sizes can be varied to accommodate different sizes of vessels.

Centre ring - removed, aperture 30mm diameter hole

Ring 1 - removed, aperture 42mm diameter hole

Ring 2 - removed, aperture 58mm diameter hole

Ring 3 - removed, aperture 78mm hole

Ring 4 - removed, aperture 103mm hole

Ring 5 - removed, aperture 128mm hole

Ring 6 - removed, aperture 153mm hole



| ACCESSORY        | 1026E  |
|------------------|--|
| Energy Regulator | ER - 1<br>Can be connected to the boiling bath, allowing control of the amount of boiling from a gentle simmer to a vigorous boil. |



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

