

1100 Series Boiling Baths

Rectangular, aluminium and stainless steel electrically heated, multi-place baths with fixed temperature for boiling water, fitted with resettable cut-out for safety and concentric ringed lids to vary the hole size.

The ER range regulates from a gentle simmer to a vigorous boil to whereas the E range has a timed temperature.



	1105E	1105ER	1101E
GENERAL			
Temp. Range	Boiling at 100°C	Simmer to boiling at 100°C	Boiling at 100°C
Temp. Control	Constant boiling no control	Energy regulator controls the amount of boiling from a gentle simmer to vigorous boil	Constant boiling no control
Heater Power (220-240V)	1000W		
Voltage	230V models 50-60Hz		
DIMENSIONS			
Tank Capacity	5 litres (Nominal Capacity)		4.5 litre
External (mm)	400w x 200d x 90h		321w x 219d x 90h

	1101ER	1104E	1104ER
GENERAL			
Temp. Range	Simmer to boiling at 100°C	Boiling at 100°C	Simmer to boiling at 100°C
Temp. Control	Energy regulator controls the amount of boiling from A gentle simmer to vigorous boil	Constant boiling no control	Energy regulator controls the amount of boiling from A gentle simmer to vigorous boil
Heater Power (220-240V)	1000W		
Voltage	230V models 50-60Hz		
DIMENSIONS			
Tank Capacity	4.5 litre	8.8 litre	
External (mm)	321w x 219d x 90h	425w x 321d x 90h	

By removing concentric rings, the hole sizes can be varied to accommodate different sizes of vessels simultaneously.

Number of holes in bath lid 6 or 12 off with maximum hole diameter 80mm.

In each aperture a 4 ring set is fitted.

Centre ring - removed, aperture 30mm diameter hole

Ring 1 - removed, aperture 42mm diameter hole

Ring 2 - removed, aperture 58mm diameter hole

Ring 3 - removed, maximum diameter 80mm hole



ACCESSORY	1105E	1101E	1104E
Energy Regulator	ER - 1 Can be connected to the boiling bath, allowing control of the amount of boiling from a gentle simmer to a vigorous boil.		



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

