

NE3DT Series Digital Water Baths

Clifton digital NE3DT series baths come complete with run back and delay start timer, clear hinged polycarbonate lid and feature PID precision temperature control.

Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance.



NE3-28DT with Polycarbonate Lid



NE3 - 4DT

NE3 - 8DT

NE3 - 14DT

GENERAL			
Temp. Range	Ambient + 5°C - 100°C		
Temp. Alarm	Over temperature alarm is automatically set 4°C above set point		
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C		
Temp. Control & Run Back Timer	Visual indicator and end message 00:01 to 99:59 (hh:mm) activates at set temperature		
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours		
Sensitivity/Uniformity	±0.2°C / ±0.1°C		
Heater Power (220-240V)	400W	800W	1000W
Drain (12mm Bore)	Drain tap fitted		
Voltage	120V or 230V models available 50-60Hz		
DIMENSIONS			
Tank Capacity	4 Litre	8 Litre	14 Litre
External (mm)	332w x 185d x 290h	332w x 270d x 290h	361w x 332d x 290h
Internal (mm)	300w x 150d x 150h	300w x 240d x 150h	325w x 300d x 150h
Water Level Max. (mm)	130		
Water Level Min. (mm)	55		

NE3 - 28DT

NE3 - 56DT

GENERAL			
Temp. Range	Ambient + 5°C - 100°C		
Temp. Alarm	Over temperature alarm is automatically set 4°C above set point		
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C		
Temp. Control & Run Back Timer	Visual indicator and end message 00:01 to 99:59 (hh:mm) activates at set temperature		
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours		
Sensitivity/Uniformity	±0.2°C / ±0.1°C		
Heater Power (220-240V)	1500W	2500W	
Drain (12mm Bore)	Drain tap fitted		
Voltage	120V or 230V models available 50-60Hz		
DIMENSIONS			
Tank Capacity	28 Litre	56 Litre	
External (mm)	538w x 332d x 290h	660w x 540d x 290h	
Internal (mm)	500w x 300d x 200h	620w x 500d x 200h	
Water Level Max. (mm)	180		
Water Level Min. (mm)	55		

Accessories on Page 14-17

Unstirred Bath Accessories - Lids

	Stainless Steel Gabled Lid with Insulated Handle	Stainless Steel Hinged Gable Lid with Insulated Handle	Stainless Steel Flat One Piece Lid with Insulated Handle	Clear Lift Off Gable Lid
UNSTIRRED BATHS				
NE1-2.5	Not Available	Not Available	LD-2.5	Not Available
NE1-2/2/2, NE1D-2/2/2, NE2-2D	SL1-4	Not Available	LD-4	GL1-4AF
NE1-4, NE1D-4, NE1-4/4/4, NE1D- 4/4/4, NE2-4D, NE2D-4/4	SL1-4	Not Available	LD-4	GL1-4AF
NE2-5D	SL1-14	Not Available	LD-14	GL1-14AF
NE1-8, NE1D-8, NE2-8D, NE2D-4/8, NE2D-8/22	SL1-8	Not Available	LD-8	GL1-8AF
NE1-9, NE2-9D	SL1-22	SL1-22H	LD-22	GL1-22AF
NE2-10D	SL1-8	Not Available	LD-8	GL1-8AF
NE1-14, NE1D-14 NE2-14D	SL1-14	SL1-14H	LD-14	GL1-14AF
NE1B-14, NE1B-18W	Not Available	SL1-14H/SS	LD-14	GL1-14AF
NE1-22, NE1D-22, NE2-22D, NE2D-8/22	SL1-22	SL1-22H	LD-22	GL1-22AF
NE1-28, NE2-28D	SL1-22	SL1-22H	LD-22	GL1-22AF
NE1-56, NE2-56D	SL1-56	Not Available	LD-56	Not Available
NE1B-28	SL1-22SS	Not Available	LD-22	GL1-22AF





<p>Stainless Steel Lid with Concentric Rings SLR1-4 : 2 x 105 mm holes SLR1-8 : 4 x 83mm holes SLR2-8 : 2 x 105 mm holes SLR1-14 : 4 x 105mm holes SLR2-14 : 6 x 83mm holes SLR3-14 : 4 x 114.5mm holes SLR1-22 : 6 x 105mm holes SLR3:22 : 6 x 114.5mm holes</p>	<p>Clear Hinged Gable Lid</p>	<p>Polypropylene Spheres: Pack of 200</p>
		
Not Available	Not Available	BP0368 : 1 Pack
SLR1-4	GL1-4	BP0368 : 1 Pack
SLR1-4	GL1-4	BP0368 : 1 Pack
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-8 and SLR2-8	GL1-8	BP0368 : 2 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 2 Packs
SLR1-8 and SLR2-8	GL1-8	BP0368 : 2 Packs
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs
Not Available	Not Available	BP0368 : 4 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs

Accessories for the NE1 B-18W Boiling Bath

ST-1	Raise and Lower Sterilising Tray
SA01409	Pair of Guide Rail Assembles for ST-1 Tray
SL1-14H/SS	Stainless Steel Hinged Gable Lid to suit 14 litre and 18 litre Boiling Baths
GL1-14AF	Lift Off Gable Lid to suit 14 and 18 litre Boiling Baths
GL1-14	Clear Hinged Gable Lid to suit 14 and 18 litre Boiling Baths

Unstirred Bath Accessories

	Stainless Steel Test Tube Racks	PVC Churchill Racks	Drain Tap
	6870 : 26 holes x 17mm DIA 6871 : 16 holes x 26mm DIA 6872 : 36 holes x 13mm DIA 6873 : 18 holes x 19mm DIA 6900 : 12 holes x 32mm DIA	C18-4 : 28 holes x 18mm DIA C18-8 : 44 holes x 18mm DIA C18-14 : 68 holes x 18mm DIA C25-4 : 28 holes x 25mm DIA C25-8 : 44 holes x 25mm DIA C25-14 : 68 holes x 25mm DIA	
UNSTIRRED BATHS			
NE1-2.5	1 Rack From Above Range	Not Available	Not Available
NE1-2/2/2, NE1D-2/2/2, NE2-2D	1 Rack From Above Range Chamber	Not Available	Not Available
NE1-4, NE1D-4, NE1-4/4/4, NE1D-4/4/4, NE2-4D, NE2D-4/4	1 Rack From Above Range	C18-4 and C25-4	Not Available
NE2-5D	4 Racks From Above Range	Not Available	Not Available
NE1-8, NE1D-8 NE2-8D, NE2D-4/8, NE2D-8/22	2 Racks From Above Range	C18-4 and C25-8	Not Available
NE1-9, NE2-9D	6 Racks From Above Range	Not Available	Not Available
NE2-10D	2 Racks From Above Range	Not Available	Not Available
NE1-14, NE1D-14 NE2-14D	4 Racks From Above Range	C18-14 and C25-14	Not Available
NE1B-14	4 Racks From Above Range	Not Available	Not Available
NE1B-18W	4 Racks From Above Range	Not Available	BX0953
NE1-22, NE1D-22, NE2-22D, NE2D-8/22	6 Racks From Above Range	Not Available	BX0953
NE1B-28, NE1-28, NE2-28D	6 Racks From Above Range	Not Available	BX0953
NE1-56, NE2-56D	12 Racks From Above Range	Not Available	BX0953

<p>Flat n Fold™ Test Tube Racks</p> <p>6770 : 26 holes x 17mm DIA 6771 : 16 holes x 26mm DIA 6772 : 36 holes x 13mm DIA 6773 : 18 holes x 19mm DIA 6800 : 14 holes x 32mm DIA</p>	<p>Stainless Steel Wire Baskets</p>	<p>Bath Armor™</p>	<p>Thermometer Clip With Bent Stem Thermometer</p>
			
1 Rack From Above Range	Not Available	BA-5 : 4 Litre Pack	TC-1
1 Rack From Above Range Chamber	Not Available	BA-5 : 4 Litre Pack	TC-1
1 Rack From Above Range	BSK-4	BA-5 : 4 Litre Pack	TC-1
4 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
2 Racks From Above Range	BSK-8	BA-12 : 12 Litre Pack	TC-1
6 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
2 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
4 Racks From Above Range	BSK-14	BA-12 : 12 Litre Pack	TC-1
4 Racks From Above Range	BSK-14	Not Available	TC-1
4 Racks From Above Range	Not Available	Not Available	TC-1
6 Racks From Above Range	BSK-22	BA-12 : 2 x 12 Litre Pack Required	TC-1
6 Racks From Above Range	Not Available	BA-12 : 2 x 12 Litre Pack Required	TC-1
12 Racks From Above Range	Not Available	BA-12 : 2 x 12 Litre Pack Required	TC-1

Common Accessories

BX0616	Drain Syphon
BKF-1	Bar Keepers Friend



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.