

Cole-Parmer®

Cole-Parmer® PB-350 Series High Capacity Portable Toploading Balances

- Five-year industry-leading warranty
- Long-lasting Internal Battery
- Rugged Housing
- Economical Price



Cole-Parmer® PB-350 Series High Capacity Portable Toploading Balances

The PB-350 series of balances offer precision and portability in a rugged housing that can stand up to demanding environments. From the classroom to the field to the manufacturing floor, PB-350 series balances are built to go anywhere that you need.

Internal or External calibration options, onboard memory storage and processing, and multiple connectivity options make the PB-350 series of balances a perfect addition to any laboratory.



Ordering Information (External Calibration)

Description	Ordering No.	Series No.	Model No.
1000g 0.01g High Capacity Portable Toploading Balance	EW-10402-58	PB-350	PB-350-1000
2000g 0.01g High Capacity Portable Toploading Balance	EW-10402-59	PB-350	PB-350-2000
6000g 0.1g High Capacity Portable Toploading Balance	EW-10402-60	PB-350	PB-350-6000
10000g 0.1g High Capacity Portable Toploading Balance	EW-10402-61	PB-350	PB-350-10000
20000g 0.1g High Capacity Portable Toploading Balance	EW-10402-62	PB-350	PB-350-20000

Technical Specification (External Calibration)

Specification	PB-350-1000	PB-350-2000	PB-350-6000	PB-350-10000	PB-350-20000
Max Capacity	1000 g	2000 g	6000 g	10000 g	20000 g
Readability	0.01 g	0.01 g	0.1 g	0.1 g	0.1 g
Calibration Type	External	External	External	External	External
Linearity	0.03 g	0.03 g	0.2 g	0.3 g	0.3 g
Repeatability	0.01 g	0.01 g	0.1 g	0.1 g	0.1 g
Stabilization Time	3 s	3 s	3 s	3 s	3 s
Weigh Pan Dimensions	195 mm x 195 mm				
Device Dimensions	333 mm x 206 mm x 97 mm				
Connectivity	2×RS232, USB-A, USB-B				

Ordering Information (Internal Calibration)

Description	Ordering No.	Series No.	Model No.
1000g 0.01g High Capacity Portable Toploading Balance	EW-10402-63	PB-350	PB-350-1000i
2000g 0.01g High Capacity Portable Toploading Balance	EW-10402-64	PB-350	PB-350-2000i
6000g 0.1g High Capacity Portable Toploading Balance	EW-10402-65	PB-350	PB-350-6000i

Technical Specification (Internal Calibration)

Specification	PB-350-1000i	PB-350-2000i	PB-350-6000i
Max Capacity	1000 g	2000 g	6000 g
Readability	0.01 g	0.01 g	0.1 g
Calibration Type	Internal	Internal	Internal
Linearity	0.03 g	0.03 g	0.2 g
Repeatability	0.01 g	0.01 g	0.1 g
Stabilization Time	3 s	3 s	3 s
Weigh Pan Dimensions	195 mm x 195 mm		
Device Dimensions	333 mm x 206 mm x 97 mm		
Connectivity	2×RS232, USB-A, USB-B		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.