

Cole-Parmer®

Cole-Parmer® Paraffin Wax Dispenser - FB-200

- Temperature control as opposed to arbitrary power control
- Improved heat up rates 45°C nominal temperature for 2.25 litres water, attained within 40 minutes
- Improved temperature stability of $\pm 0.5^{\circ}\text{C}$
- Overheating prevented with new temperature control mechanism
- PTFE black interior and powder-coated white aluminium exterior
- Virtually maintenance-free and extremely easy to clean
- Black interior provides high contrast for viewing samples



Cole-Parmer® Paraffin Wax Dispenser - FB-200



FB-200 Paraffin Wax Dispenser

Cole-Parmer® paraffin section flotation bath is designed to assist with the handling of paraffin wax samples in histology and pathology laboratories. It is essentially a hot distilled water floating out bath that allows for the meticulous manipulation and location of sections onto glass slides.

The paraffin section flotation bath has improved heat up times achieving a stable 45°C using 2.25 litres of water within 40 minutes. The stability and accuracy have also been improved due to the change of control which now also allows the user to set the temperature as opposed to the arbitrary power, thereby, reducing the risk of overheating sections.

The paraffin section flotation bath has a PTFE black interior and powder-coated aluminium exterior, which makes it virtually maintenance-free, extremely easy to clean and allows for improved visibility of floating sections.

Technical Specification

Specification	FB-200
Capacity	2.5L
Temperature	30°C to 70°C
Temperature Stability	±0.5°C
Bath material	Aluminium base with PTFE coated interior
Dimensions (d x w x h)	330 x 330 x 140 mm
Weight	2.7kg
Dimensions, mm (w x d x h)	240W

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Paraffin Section Flotation Bath 230V UK & EU	48572-55	FB-200	FB-200	MH8517
Paraffin Section Flotation Bath 115V US	48572-50	FB-200	FB-200-115	MH8517X1
Paraffin Bath Lid	48572-99	-	-	AZ9241-



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.