


Cole-Parmer®

Cole-Parmer® Tube Rollers, RS-200 and RS-200D

- Greater capacity and multiple rollers saves bench top space and meets standard lab throughput needs
- Featuring removable side panels
- Featuring detachable roller tubes enabling large volume/diameter bottles to fit
- Choice of analog or digital control
- Choice of fixed or variable speed



6 Tube Roller RS-200

Cole-Parmer® Tube Rollers, RS-200 and RS-200D

Cole-Parmer® roller mixers provide a gentle but efficient mixing ideal for mixing blood samples, viscous substances and liquid-solid suspensions where minimum aeration is required. Cole-Parmer® is committed to providing the highest levels of usability and as such has added additional features to the roller mixers. RS-200 roller mixers now feature easy to remove side panels. The side panels provide the user with complete peace of mind that their sample will be mixed effectively during continued use. The side panels ensure the sample bottles do not slide off the mixer and therefore provide the highest levels of both safety and functionality.

The Cole-Parmer® roller tubes provide the highest levels of versatility allowing the roller mixers to conveniently detach to accommodate larger vessels and therefore allow the user to mix a wide variety of differently sized tubes and bottles.



6 Tube Roller RS-200D

Technical Specification

Specification	RS-200	RS-200D
Number of rollers	6/9	6/9
Speed	33rpm	5 to 60rpm
Amplitude	16mm	16mm
Maximum load	3.0kg	3.0kg
Controls	Analog	Digital
Timer	No	Yes
Roller size (length x diameter)	340x30mm	340x30mm
Overall dimensions (w x d x h)	595x255x110mm 565x360x110mm	595x255x110mm 565x360x110mm
Net weight	5.1kg & 6.9kg	5.1kg & 6.9kg
Electrical Supply	230V, 50Hz, 50W 120V, 60Hz, 50W	230V, 50Hz, 50W 120V, 60Hz, 50W

Ordering Information

Description	Ordering Number	Series Number	Model Number	Legacy Cole-Parmer® SKU	Legacy Stuart® SKU
Roller mixer, 6 rollers, analog control, fixed speed, 230V 50 Hz	51901-18	RS-200	RS-200	04653-02	SRT6
Roller mixer, 6 rollers, analog control, fixed speed, 120V, 60Hz	51901-19	RS-200	RS-200-120	04653-01	SRT6/120V/60
Roller mixer, 6 rollers, digital control, variable speed, 230V 50 Hz	51901-22	RS-200	RS-200D	04653-06	SRT6D
Roller mixer, 6 rollers, digital control, variable speed, 120V, 60Hz	51901-23	RS-200	RS-200D-120	04653-05	SRT6D/120V/60
Roller mixer, 9 rollers, analog control, fixed speed, 230V 50 Hz	51901-20	RS-200	RS-200-L	04653-04	SRT9
Roller mixer, 9 rollers, analog control, fixed speed, 120V, 60Hz	51901-21	RS-200	RS-200-L-120	04653-03	SRT9/120V/60
Roller mixer, 9 rollers, digital control, variable speed, 230V 50 Hz	51901-24	RS-200	RS-200D-L	04653-08	SRT9D
Roller mixer, 9 rollers, digital control, variable speed, 120V, 60Hz	51901-25	RS-200	RS-200D-L-120	04653-07	SRT9D/120V/60



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.