

# Cole-Parmer®

## Cole-Parmer® Stirrer Multiposition, ST-200-MP

- Three independently controlled stirring positions
- Powerful magnets for strong magnetic coupling
- Stainless steel top plate





## Cole-Parmer® Stirrer Multiposition, ST-200-MP

A space saving and economical unit with three independent stirring positions in a footprint only 600 x 270mm. Easily accommodates 3 x 2 litre beakers. The stainless steel top plate is robust and, unlike aluminium, does not produce eddy currents and so ensures a very powerful coupling and stirring action. Powerful magnets and motor give stirring speed up to 1500rpm and volumes up to 15 litres \*. Supplied complete with 3 x 25mm PTFE coated stir bars.

### Technical Specification

Specification	ST-200-MP
Plate material	Stainless steel
Plate dimensions, mm	550 x 210
Speed range, rpm	100 to 1500
Dimensions, mm (w x d x h)	606 x 268 x 94
Net weight, kg	7
Electrical supply	230V, 50Hz, 150W
IP Rating	31

### Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Stirrer, magnetic, 3 position and Allen key	51400-02	ST-200	ST-200-MP	SB161-3
Retort rod bracket	04805-96	ST-200	N/A	SR3
Retort rod, 600mm x 12mm diameter.	04805-91	ST-200	N/A	SR1



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.