



Dutscher Group
 Sonnenschein 8
 D-48565 Steinfurt
 Tel: (+49) 02551-864310
 Fax: (+49) 02551-864312
 contact@kisker-biotech.com
 www.kisker-biotech.com

Specification

Product name: H-ATDA - Adhesive crepe tapes

REF	Width x Length
H-ATDA	19 mm x 50 m

Materials and Properties

- Backing:** Crepe paper
- Indicator:** Steam indicator according to EN ISO 11140-1, class 1.
 Colour change from white to dark grey – black.
- Colour:** beige
- Adhesive:** normal adhesive
- Total thickness:** 0,165 mm
- Adhesion to steel:** 8 N / 25 mm
- Tensile strength:** 90 N / 25 mm
- Elongation at break:** 7 %
- Toxicity:** toxic effects are not known

Compatibility:

- Packaging material:**
- | | | |
|-------------------|--|--|
| Crepe Paper | <input checked="" type="checkbox"/> suitable | <input type="checkbox"/> unsuitable |
| Fleece | <input checked="" type="checkbox"/> suitable | <input type="checkbox"/> unsuitable |
| Textiles material | <input type="checkbox"/> suitable | <input checked="" type="checkbox"/> unsuitable |
- Sterilization methods:**
- | | | |
|---------------------|--|--|
| STEAM | <input checked="" type="checkbox"/> suitable | <input type="checkbox"/> unsuitable |
| EO (Ethylene oxide) | <input type="checkbox"/> suitable | <input checked="" type="checkbox"/> unsuitable |
| FORM (Formaldehyde) | <input type="checkbox"/> suitable | <input checked="" type="checkbox"/> unsuitable |
| VH2O2 (Plasma) | <input type="checkbox"/> suitable | <input checked="" type="checkbox"/> unsuitable |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

