

REFERENCE :
062583CL

ClearLine micofuge for PCR plates

Unique rotor design prevents sample spillage

Plate support at 75° before and after centrifugation


90° plate support during centrifugation

Fixed speed

Centrifugation stops when the lid is opened

Shutdown in 4 seconds

Compatible with all types of PCR plates (with skirt, semi or not), microplates as well as Deepwell plates (with optional rotor)



Before/After centrifugation During centrifugation

Speed	2550 rpm / 600g
Capacity	2 PCR plates* or 2 microplates or 192 tubes PCR 0.2ml or 24 strips of 8 tubes
Dimensions L x l x H (mm)	230 x 260 x 197
Weight (kg)	4.1
Electrical	230V, 50-60 Hz
Warranty	2 years

*Maximum capacity of 4 PCR plates when stacked with the Deepwell plate rotator

Reference	Description
062584	Adapter for 0.2ml PCR tubes (x 2)
062585	Rotor for Deepwell plates (max. height 35mm)

Units :

1 centrifuge



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.