

## Clearline Micropipette

### Features :

- ✓ Volume locking prevents accidental counter shift
  - ✓ Ergonomic : contoured handle to adapt to the hand, finger support to reduce wrist fatigue
  - ✓ Ejector with infinitely variable length adjustment makes it easier for different tips
  - ✓ Low pipette weight
  - ✓ Soft springs system provides comfort at work
  - ✓ Fully autoclavable
  - ✓ UV resistant
  - ✓ Calibration key included
  - ✓ Each instrument passes stringent validation and control procedures in accordance with EN ISO 8655
- 
- ✓ Volume : 0,5 - 10  $\mu$ l
  - ✓ Color : red
  - ✓ Channel : singlechannel



**Units : 1 micropipette**



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.