

REFERENCE :

062650CL

**Programmable rotator Mini
ClearLine®**

Features :

- Rotating shaker for tubes from 1.5 to 50 ml
- Wide range of rotation and agitation modes
- 9 pre-recorded programs, 2 customizable programs
- 360° rotation with the possibility of pausing and shaking
- Use in cold room and incubator
- Optional platform for rocker function
- Supplied with: 2 racks for 6 x 1.5 ml tubes
2 racks for 1 tube of 15 ml
2 racks for 1 tube of 50 ml



Speed (rpm)	1 to 80
Angle	Up to 360°
Timer	1 sec to 99 h 59 min
Capacity	6 tubes of 50 ml 10 tubes of 15 ml 60 tubes of 1,5 ml
Dimensions L x P x H (mm)	290 x 137 x 199
Temperature	2 to 70°C
Weight (kg)	1,8
Power supply	230V, 50-60Hz

F1	F2	F3	F4	F5	F6	F7	F8	F9	--	C1 C2	PP*	Key:
										Custom Mode	Platform Rocking Mode	 Pause

Unité de vente :

1 agitateur



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.