

REFERENCE :

060956CL

ClearLine® Orbital shaker mini

Features :

- Orbital shaker with digital screen
- Compact
- Capacity : 4 x 1l
- Compatible with platforms « Magic Clamp » (magnetic)
- Comes with a classic platform and a slip mat 240 x 200 mm
- 2 years warranty



Speed (rpm)	50 to 300
Increment (rpm)	10
Orbit (mm)	19
Maximum load (kg)	4
Platform dimensions (mm)	240 x 200
Shaker Dimensions L x D x H (mm)	300 x 300 x 120
Display	Digital
Operating Temperature	4 to 65°C
Weight (kg)	5,5
Electrical power	230V, 50-60Hz

Sales Unit :

1 Shaker



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.