

**REFERENCE :**

**062819CL**

## ClearLine® Orbital shaker

### Features :

- Orbital shaker with touch screen
- Max capacity : 2 x 2l
- Compatible with platforms « Magic Clamp »
- Comes with a classic platform and a slip mat 330 x 300 mm



Speed (rpm)	30 to 300
Increment (rpm)	1
Orbit (mm)	19
Maximum load (kg)	6,4
Platform dimensions (mm)	330 x 300
Shaker Dimensions L x D x H (mm)	280 x 160 x 75
Display	Digital
Operating Temperature	4 to 65°C
Weight (kg)	11,25
Electrical power	230V, 50-60Hz

### Sales Unit :

1 Shaker



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.