

## Tube Roller RM 2



### Description

Gentle tilting and rotating motion for uniform mixing

#### Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

#### Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

### Details

<b>Motion:</b>	rotating and tilting motion
<b>Shaking platform:</b>	200 x 300 mm
<b>Max. load:</b>	2 kg
<b>Shaking speed:</b>	5 - 80 rpm, in steps of 5
<b>Angle of tilt:</b>	3° fixed
<b>Stroke at end of rolls:</b>	15 mm

---

**Dimensions of rolls:** 300 x Ø 32 mm  
**Number of rolls:** 6 pieces  
**Runtime:** 0 - 99 h 59 min / continuous  
**Order No.:** 6138 000  
**Accessory:** Overhead Rack RM 2

## Technical specifications

**Electrical supply:** 100 - 240 V, 50 / 60 Hz  
**Enclosure protection:** IP 21  
**Heat emission:** max. 5 W  
**Ambient temperature:** 5°C to 45°C  
**Relative humidity:** ~ 85 % (non-condensing)  
**Dimensions (w x d x h):** 244 x 465 x 115 mm  
**Weight:** 7,3 kg

---



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.