


Miniature Shaker KM 2



Description

The orbital shaker KM 2 has a shaking platform of 200 x 295 mm which can be loaded with 2 kg. The stroke is 20 mm and the shaking speed between 5 and 220 rpm.

Basic Equipment

Basic device incl. basic platform and rubber mat

Loading Capacity

Erlenmeyer flasks 100 ml 9 pieces

Erlenmeyer flasks 500 ml 2 pieces

You will find further loading possibilities [here](#)

Details

Motion:	orbital
Shaking platform:	200 x 295 mm
Max. load:	2 kg
Shaking speed:	5 - 220 rpm, in steps of 5
Stroke:	20 mm
LED display:	0 - 99 h 59 min / continuous
Order No.:	6113 000
Accessory:	Universal tray KM Mini Universal tray KM Clamping strip KM

Technical specifications

Electrical supply:	100 - 240 V, 50 / 60 Hz
Enclosure protection:	IP 21
Heat emission:	approx. 4 W
Ambient temperature:	4°C to 45°C
Relative humidity:	~ 95 % (non-condensing)
Dimensions (w x d x h):	250 x 415 x 145 mm
Weight:	8,1 kg



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.