

# Microplate Shaker TiMix 2



## Description

The Microplate Shaker TiMix 2 especially for small samples in microplates, deepwell-plates or other plates in microplate format.

### Basic Equipment

Basic device incl. basic platform with holders for 4 microplates

### Loading Capacity

4 microplates  
(deepwell-plates, or other plates in microplate format)

## Details

<b>Motion:</b>	exact orbital motion
<b>Shaking platform:</b>	200 x 295 mm for 4 microplates
<b>Max. load:</b>	2 kg
<b>Shaking speed:</b>	100 - 1400 min <sup>-1</sup> , steplessly adjustable
<b>Stroke:</b>	3 mm
<b>Timer:</b>	0 - 120 min / continuous
<b>Order No.:</b>	6110 000
<b>Accessory:</b>	Shaking platform KM Clamping strip KM

## Technical specifications

<b>Electrical supply:</b>	230 V or 115 V, 50 / 60 Hz please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	approx. 15 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	220 x 410 x 110 mm
<b>Weight:</b>	8,5kg

---



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.