

## SIEVING MACHINE AIR SIZER 200



**The new Air Sizer 200 is ideal for sieving very fine dry particles, which require efficient dispersion and desagglomeration via air jet technology (e.g. electrostatic material).**

It is also the perfect instrument to quickly provide a sieve cut of powdered materials.



[Click to view video](#)

Product Video: AIR SIZER 200

## ADVANTAGES

- | Advanced air jet technology for fine particles, usable for dry material 20 µm upwards
- | Adjustable nozzle speed, 5 - 55 rpm
- | Extremely efficient & fast sieving times
- | Sieving action keeps apertures clear
- | Air flow fluidises and helps to separate sample
- | Ideal for electrostatic materials
- | Pre-filter unit & industrial vacuum available as accessories
- | Maintenance free

## TECHNICAL DATA

<b>Range</b>	20 µm - 4 mm
<b>Drive / sieving motion</b>	dispersion by air jet
<b>Max. number of sieves</b>	1
<b>Max. mass of sieve stack</b>	a premium air jet sieve
<b>Speed</b>	5 - 55 rpm (nozzle speed)
<b>Time display</b>	digital, 0:10-99:50 min
<b>Vacuum</b>	20 - 99 mbar
<b>Suitable for dry sieving</b>	yes
<b>Suitable for wet sieving</b>	-
<b>Max. height of sieve stack</b>	1 sieve
<b>Accessories</b>	pre-filter unit / industrial vacuum
<b>Model</b>	benchtop
<b>Protection code</b>	IP 40
<b>Electrical supply</b>	100 - 240 V , 50/60Hz
<b>Power connection</b>	1 - phase
<b>W x H x D</b>	450 x 235 x 435 mm
<b>Net weight</b>	~ 16 kg
<b>Standards</b>	CE



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

