



Product Code: FSL400G

Slimline 400 Litre Upright Glass Door Refrigerated Cabinet



Features & Benefits

- When space is a prime concern, Slimline Cabinets provide the ideal storage solution
- LED electronic temperature display for accurate control and clear visual reference
- Glass door and light (25°C ambient ISO Climate Class 3)
- Bottom mounted refrigeration system draws air from the cooler part of the kitchen ensuring efficient operation and lower running costs

Shelving

Nylon Coated Shelf Option

Refrigeration

R134a - Industry Standard Standard

R290 - Energy Efficient Standard

Integral Refrigeration System Standard

Finish

s/steel ext, glass door & alu intStandard

Base Fittings

Castors 80mm Standard

Electrical & Service Data

230/50/1 Standard

Fuse Required (A) 10

Doors

Glass Standard

Left Hand Hinged Door Option

Specification

Refrigerator +1/+4°C Standard

Width (mm) 600

Depth (mm) 705

Height (mm) 1990

Gross Volume (Ltr) 400

Noise dBA @ 1m TBC

Footnotes

- Bottom mounted system draws air from the cooler part of the kitchen ensuring lower running costs
- Forced air refrigeration for efficient operation, even temperature distribution and fast recovery
- Automatic Defrost for hassle-free operation
- Electronic temperature display for accurate control
- Delivers the reliability and performance you would expect from a Foster at a highly competitive price



Energy



Quality



Value



Design



Hygiene



After Care

All images shown are for illustration purpose only. All images and photographs are not intended to be relied upon for, nor to form part of, any contract unless specifically incorporated in writing into the contract. Solely the dimensions and specifications incorporated in the quotation or purchase order confirmation shall be contractually binding.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.