



# Product Code: XR415H 33-111

## xtra by Foster slimline refrigerator

---



## Features & Benefits

---

- The xtra range: Built for you, designed for your budget
- Reliable and robust refrigeration at an affordable price
- In stock and available now for immediate shipping
- Robust stainless steel finish
- Easy read microprocessor controller provides simple and accurate temperature control
- Full Foster 2 years parts and labour warranty

## Doors

<b>Solid</b>	Standard
<b>Left Hand Hinged Door</b>	Option

## Specification

<b>Refrigerator +3/+5°C</b>	Standard
<b>Width (mm)</b>	600
<b>Depth (mm)</b>	650
<b>Height (mm)</b>	1875
<b>Gross Volume (Ltr)</b>	410

## Refrigeration

<b>R290 - Energy Efficient</b>	Standard
--------------------------------	----------

## Finish

<b>s/s ext &amp; alu int</b>	Standard
------------------------------	----------

## Electrical & Service Data

<b>230/50/1</b>	Standard
<b>220/60/1 Electrical Supply</b>	Option
<b>Fuse Required (A)</b>	10

## Energy

<b>AEC</b>	873.31
<b>Nett Volume</b>	284.8
<b>Energy Grade</b>	E
<b>Energy Efficiency Index</b>	81.09

## Footnotes

- Extra shelf and pair of trayslides available as an option
-

# Energy Performance Data and information



**Foster Refrigerator** - a division of ITW Ltd  
Oldmedow Road, KING'S LYNN, Norfolk PE30 4JU  
Tel: 0843 216 8800 Fax: 0843 216 4700  
Registered in England No. 00559693  
VAT No. GB 792 4414 16



Master data version: Ver.01.0014AB

Generated: 31/05/2019

## Application

This appliance is intended for use in ambient temperatures up to {1}°C

## Safety

Precautionary instructions and maintenance advice, to ensure efficient operation, can be found in this products operating instructions.

## Recycling

For recycling and end of life care for this product, please refer to environmental section of the website or visit [our policies](#)



Energy



Quality



Value



Design



Hygiene



After Care

All images shown are for illustration purpose only. All images and photographs are not intended to be relied upon for, nor to form part of, any contract unless specifically incorporated in writing into the contract. Solely the dimensions and specifications incorporated in the quotation or purchase order confirmation shall be contractually binding.



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.