


| | | | | | |
|--------------------------------|--------------------------------------|-------------------------|-------------------------|-----------------------------|--|
| PDS No. 6071xx / 607160-TRI | PRODUCT DATA SHEET | | | | Page 1 of 1 |
| Revision 09 | Pipettes, CELLSTAR®, 10 ml | | | |  greiner bio-one |
| | Greiner Item-No. 6071xx / 607160-TRI | | | | |
| Valid for Item-No.: | 607107 (sterile) | 607160 (sterile) | 607180 (sterile) | 607160-TRI (sterile) | |

| 1. | Description / Specification | |
|-----|-----------------------------|---|
| 1.1 | Description | PS Pipette, graduation (scale 1/10 ml), lateral orange stripe, with synthetic plug, 10 ml, sterile |
| 1.2 | Dimensions | See Customer Drawing |
| 1.3 | Volume | 10 ml; volumetric accuracy $\pm 2\%$ |
| 1.4 | Material / Resin | Pipette and nozzle: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal |
| 1.5 | Colour | Pipette: clear Plug: white |
| 1.6 | Sterilization | 607107, -160, -180: SAL 10^{-3} 607160-TRI: SAL 10^{-6} |
| 1.7 | Quality Control | - <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information | For single use only |

| 2. | Features | |
|-----|--------------------------|--|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic. |
| 2.2 | Temperature range | 0°C to +40°C |
| 2.3 | Autoclavability | No |
| 2.4 | Centrifugation, max. RCF | N/A |
| 2.5 | Chemical Resistance | See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html |
| 2.6 | Shelf life | 5 years after month of production |
| 2.7 | Other Information | - |

| 3. | Packaging | 607107 | 607160 | 607180 | 607160-TRI |
|-----|-------------------|--|--|--|--|
| 3.1 | Pieces / Bag | 25 Reclosable bulk packs (plastic) | 1 "peel off"-packaging (plastic/plastic) | 1 "peel off"-packaging (paper/plastic) | 1 "peel off"-packaging (plastic/plastic) 10 "Triple packed" |
| 3.2 | Pieces / Box | 500 | 200 | 200 | 100 |
| 3.3 | Lot-No. | F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) | F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) | F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality | | | |

| 4. | Other Information |
|-----|-------------------|
| 4.1 | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|-------------------------|-------------------------|-------------------------|---|
| Revision 08 | Date 30 January 2019 | Date 31 January 2019 | Date 31 January 2019 | |
| Date 01.10.2015 | Name S. Kaelberer | Name Dr. J. Bischoff | Name A. Illig | |

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

