


PDS No. 7601xx / 760160-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 09	Pipettes, CELLSTAR®, 25 ml				 greiner bio-one
	Greiner Item-No. 7601xx / 760160-TRI				
Valid for Item-No.:	760107 (sterile)	760160 (sterile)	760180 (sterile)	760160-TRI (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 2/10 ml), lateral red stripe, synthetic plug, 25 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	25 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	Pipette and nozzle: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Plug: white
1.6	Sterilization	760107, -160, -180: SAL 10 ⁻³ 760160-TRI: SAL 10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	760107	760160	760180	760160-TRI
3.1	Pieces / Bag	25 Reclosable bulk packs (plastic)	1 "peel off"-packaging (plastic/plastic)	1 "peel off"-packaging (paper/plastic)	1 "peel off"-packaging (plastic/plastic) 10 "Triple packed"
3.2	Pieces / Box	200	200	200	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality			

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 08	Date 30 January 2019	Date 31 January 2019	Date 31 January 2019	
Date 09.11.2017	Name S. Kaelberer	Name Dr. J. Bischoff	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

