

Haier

UFR106 Infrared Thermometer Introduction

额温枪
国内定制专用校准文件

UFR106-- Product Features

Haier

UFR106



Measure
In 1s

Automatic
power-off

20
memories

Switch
between
Celsius and
Fahrenheit



Measure in 1 second

—
No contact with human skin, one second temperature measurement



Green Backlight

In low light, it can also clearly display the screen information, so that you can easily read the data.

UFR106--Product Parameters

Haier




ITEM	UFR106
Measurement method	Non-contact
Measurement range	32.0°C~43.0°C (89.6°F~109.4°F)
Accuracy	±0.2°C/±0.4°F (within 36°C~39°C/96.8°F~102.2°F)
Display resolution	0.1°C/0.1°F
Memory function	20 sets memory of measurement values
Power source	2pcs AAA alkaline battery Notice: no battery
Main unit weight	Approx. 75g (batteries not included)
Main unit size	L151mm * W53mm * H41mm





ITEM	UFR106
BOX SIZE	93*53*176mm
CTN SIZE	392*290*376mm
CTN	40pcs
G.W	6.5KG
N.W	5.9KG

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFIKAT ◆ 認證證書 ◆ CERTIFICATE ◆ CERTIFICATE ◆ CERTIFICATE



General Conformity Certificate for
 Compliance with Council
 Directive 90/269/EEC on the
 Approximation of the Laws,
 Regulations, Administrative
 Provisions and
 Directives of the Member States
 ZUG-95-244.10.08

Product Service

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFIKAT ◆ 認證證書 ◆ CERTIFICATE ◆ CERTIFICATE ◆ CERTIFICATE

EC Certificate

Full Quality Assurance System
 Directive 90/269/EEC on Medical Devices (MDD), Annex II excluding (4)
 (Devices in Class IIa, IIb or III)
No. G1 078672 0014 Rev. 00

Manufacturer: [Redacted]


Product Category(ies): Digital Blood Pressure Monitors, Infrared Thermometer, Compressor Nebulizers

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of the Directive and is subject to periodical surveillance. For marking of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.: G21917901

Valid from: 2020-02-21
Valid until: 2022-02-07

Date: 2020-02-21


 Christoph Dicks
 Head of Certification/Notified Body

Page 1 of 2
 TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123
 TÜV SÜD Product Service GmbH • Certification Body • Riederstraße 65 • 80399 Munich • Germany

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFIKAT ◆ 認證證書 ◆ CERTIFICATE ◆ CERTIFICATE ◆ CERTIFICATE



General Conformity Certificate for
 Compliance with Council
 Directive 90/269/EEC on the
 Approximation of the Laws,
 Regulations, Administrative
 Provisions and
 Directives of the Member States
 ZUG-95-244.10.08




Product Service

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFIKAT ◆ 認證證書 ◆ CERTIFICATE ◆ CERTIFICATE ◆ CERTIFICATE

EC Certificate

Full Quality Assurance System
 Directive 90/269/EEC on Medical Devices (MDD), Annex II excluding (4)
 (Devices in Class IIa, IIb or III)
No. G1 078672 0014 Rev. 00

Facility(ies): [Redacted]


Page 2 of 2
 TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123
 TÜV SÜD Product Service GmbH • Certification Body • Riederstraße 65 • 80399 Munich • Germany

中华人民共和国医疗器械注册证

注册证编号: [Redacted]

注册人名称	[Redacted]
注册人住所	[Redacted]
生产地址	[Redacted]
产品名称	红外体温计
型号、规格	UFR101, UFR102, UFR103, UFR105, UFR106
结构及组成	由外壳、红外传感器、微电脑控制电路及液晶显示屏组成。
适用范围	通过测量额头的热辐射来显示被测对象的体温。
附件	产品技术要求。
其他内容	无
备注	

审批部门: 广东省食品药品监督管理局 批准日期: 2018年09月08日
 有效期至: 2021年09月08日





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.