

Haier Biomedical

Intelligent Protection of Life Science

Combined Refrigerator and Freezer

HYCD-282A



Refrigerator & Freezer Combination

Independent display for refrigerator section and freezer section



Alarm Functions

High/low temperature, power failure, sensor error, low battery, high/low ambient temperature, and remote alarm



Standard Stainless Steel Interior

Both refrigerator and freezer chamber are stainless steel interior, durable, easy to clean and anti-corrosion

Ergonomic Design



Lockable design



Intelligent alarm



Foot/Caster design



Test hole design

Specifications

HYCD-282A

Technical Data

| | |
|---------------------|---|
| Type | Combined Refrigerator and Freezer |
| Cabinet Type | Upright |
| Climate Class | 4 |
| Cooling Type | Refrigerator: Forced air cooling Freezer: Direct cooling |
| Defrost Mode | Refrigerator: Auto Freezer: Manual |
| Refrigerant | CFC-Free |
| Noise Level (dB(A)) | 43 |

Performance

| | |
|-----------------------|---------------------------------------|
| Temperature Range(°C) | Refrigerator: 2~8 Freezer: -20~-40 |
|-----------------------|---------------------------------------|

Control

| | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

Electrical Data

| | |
|-----------------------|------------|
| Power Supply(V/Hz) | 220~240/50 |
| Power(W) | 400 |
| Electrical Current(A) | 2.8 |

Dimensions

| | | |
|-------------------------------|---|--|
| Capacity(L/Cu.Ft) | Refrigerator:185/6.5 | Freezer:97/3.4 |
| Net/Gross Weight(approx) | 145/160(kg) | 319.7/352.7(lbs) |
| Interior Dimensions(W*D*H) | Refrigerator:605*510*720(mm) Freezer:515*465*440(mm) | 23.8*20.0*28.3(in) 20.3*18.3*17.3(in) |
| Exterior Dimensions(W*D*H) | 736*660*1810(mm) | 29.0*26.0*71.3(in) |
| Packing Dimensions(W*D*H) | 775*720*1930(mm) | 30.5*28.3*76(in) |
| Container load (20'/40'/40'H) | 21/45/45 | |



Functions

| | |
|----------------------|---|
| High/Low Temperature | Y |
| Remote Alarm | Y |
| Power Failure | Y |
| Sensor Error | Y |
| Low Battery | Y |
| High Ambient Temp | Y |

Accessories

| | |
|-----------------|-----------------------------------|
| Caster | Y |
| Foot | Y |
| Test Hole | Y |
| Shelves/Drawers | Refrigerator: 3/1 Freezer: 0/2 |
| USB Interface | Y |

Others

| | |
|-------------|----|
| Certificate | CE |
|-------------|----|



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.