

## Pharmacy Refrigerators HYC-118



**Audible and  
Visual Alarms**  
Audible buzzer and visible  
flashing light



**Five Alarm  
Functions**  
High/low temperature,  
power failure, door ajar,  
sensor error, low battery



**Air Cooling,  
Auto Defrost**  
Better cabinet temperature  
uniformity, auto defrost

### Ergonomic Design



High Accuracy  
Temperature Control



Interior  
LED



Lockable  
Design

# Specifications

## HYC-118

### Technical Data

Cabinet Type	Upright, solid door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	32

### Performance

Temp Range(°C)	2~8
----------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	135	75
Electrical Current(A)	1.03	1.11

### Dimensions

Capacity(L/Cu.Ft)	118/4.2
Net/Gross Weight(approx)	41/46(kg) 90.4/101.4(lbs)
Interior Dimensions(W*D*H)	515*415*630(mm) 20.3*16.3*24.8(in)
Exterior Dimensions(W*D*H)	597*635*835(mm) 23.5*25.0*32.9(in)
Packing Dimensions(W*D*H)	680*690*910(mm) 26.8*27.2*35.8(in)
Container Load(20'/40'/40'H)	48/102/102



### Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

### Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	3/1

### Others

Certificate	CE,UL
-------------	-------



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.