

Pharmacy Refrigerators HYC-1378



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Special Designed Slot for Sign Card

Adjustable shelves, the slot design be clear,
store items is easy to distinguish



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



Self Closing Door

Energy conserving , maintains temperature
uniformity to protect valuable samples

Specifications

HYC-1378

Technical Data

Cabinet Type	Upright, double door		
Climate Class	N		
Cooling Type	Forced air cooling		
Defrost Mode	Auto		
Refrigerant	CFC-Free		
Noise((dB(A))	45	45	48

Performance

Temp Range(°C)	2~8		
------------------	-----	--	--

Control

Controller	Microprocessor		
Display	LED		

Electrical Data

Power Supply(V/Hz)	220-240/50	208~230/60	115/60
Power(W)	800	800	950
Electrical Current(A)	4.2	4.2	11

Dimensions

Capacity(L/Cu.Ft)	1378/48.7		
Net/Gross Weight(approx)	310/360(kg) 683.4/793.7(lbs)		
Interior Dimensions(W*D*H)	1320*700*1500(mm) 52.0*27.6*59.1(in)		
Exterior Dimensions(W*D*H)	1440*925*1980(mm) 56.7*36.4*78.0(in)		
Packing Dimensions(W*D*H)	1570*1010*2120(mm) 61.8*39.8*83.5(in)		
Container Load(20'/40'/40'H)	6/12/12		



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	12/-
USB Interface	Y
Temp Recorder	Optional

Others

Certificate	CE
-------------	----



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.