

Blood Bank Refrigerators

HXC-158B



**Infinite Divided
Stainless Steel Drawer**
Stainless steel drawer can be unlimited separated based on the demand



**Three Types
of Alarm**
Audible buzzer, visible flashing light, remote alarm



Five Alarm Functions
High/low temperature, power failure, sensor error, low battery, door ajar

Ergonomic Design



Famous
Compressor



Stainless
Steel Drawer



Condensation
Free Design



Intelligent Alarm

Specifications

HXC-158B

Technical Data

Type	Drawer-Type
Climate Class	ST
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	42

Performance

Temperature Range(℃)	2~6
Ambient Temperature(℃)	10~38

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/HZ)	220~240/50/60
Power(W)	350W
Electrical Current(A)	2.6

Dimensions

Capacity(L/Cu.Ft)	158/5.6
Blood Storage Capacity (450ml blood bags)	84
Net/Gross Weight(approx)	113/126(kg) 249.1/277.8(lbs)
Interior Dimension(W*D*H)	460*370*950(mm) 18.1*14.6*37.4(in)
Exterior Dimension(W*D*H)	560*570*1530(mm) 22.0*22.4*60.2(in)
Packing Dimension(W*D*H)	645*675*1680(mm) 25.4*26.6*66.1(in)
Container Load(20'/40'/40'H)	27/54/54



Functions

High/Low Temperature	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote Alarm	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	-
Drawers/Inner doors	4/-
USB Interface	Optional
Temp Recorder	Y

Others

Certificate	CE
-------------	----



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.