



## Arctic Squares® Polycarbonate Cryostorage Boxes

### Safely store cryogenic vials

Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease.

- Alphanumeric labeling on lid or base (depending on the unit) for quick location and retrieval of tubes
- Magnet on top of forceps enables easy storing on metal surfaces
- Stackable, saves valuable space in the lab
- Liquid nitrogen vent and drainage holes



US Patent No.: D529,622

EP Registered Design No.: 00551429-0001, 00551429-0002, 00551429-0003

Item No.	Color	Tube Array	Tubes mL	Wells	Well Dia. mm	Lid	Packaging	L x W x H in	L x W x H cm	Forceps
HS23325A	Ruby Red	5 x 5	2.0	25	12.5	Slide On	8	3 x 3 x 2	7.6 x 7.6 x 5.3	
HS23381B	Assorted*	9 x 9	2.0	81	12.5	Hinged	4	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓
HS23300D	Arctic Blue	10 x 10	2.0	100	12.5	Hinged	4	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓

\*Assorted colors (Ruby Red, Amethyst Purple, Arctic Blue, Emerald Green)



## Arctic Squares® Upright Freezer Rack

### Innovative horizontal drawer system for upright freezers

- Constructed of corrosion-resistant stainless steel with four drawers to hold 16 boxes
- Drawers feature a unique stop, that when fully extended pivot down allowing for easy box access
- Steel drawer handles provide a user-friendly grip
- Freezer rack also includes a pull-out wire handle for easy removal off the freezer shelf
- Compatible with most standard 2 in (5.3 cm) freezer storage boxes

Unique drawer stop that pivots down allowing for easy box access

Steel drawer handles have a space to insert a label, and provide a user-friendly grip

Item No.	Config.	Box Height in	Box Height cm	L x W x H in	L x W x H cm
HS2864	4 x 4	2	5.3	5.4 x 22.4 x 9.3	13.9 x 57.1 x 23.8



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.