




Cool Containers

Durable containers for chilling temperature-sensitive samples

Suitable for low temperature use. Ideal for ice, ice slurries and CO₂ pellets

- Made from polyurethane with metal inserts for stability
- Close-fit lid keeps evaporation to a minimum and features a comfortable recessed knob
- Tapered design for easy stacking
- Precise pour spout and large handles



Lid nests under bucket
for convenience and to save space

Item No.	Color	Style	Capacity liters	L x W x H in	L x W x H cm
HS28714K	Black	Bucket	4	7.1 (H) x 9.2 (D)	18.1 (H) x 23.4 (D)
HS28721G	Green	Pan	1	10.5 x 7.7 x 4.5	26.7 x 19.6 x 11.4
HS28724B	Blue	Pan	4	15 x 10 x 5.5	38.1 x 25.4 x 14
HS28729R	Red	Pan	9	20.5 x 11.5 x 6.4	52.1 x 29.2 x 16.3



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.