



# XT Series

**Extended Time Cryogenic  
Small Freezers**

# 3XTL

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS                | 3XTL         |
|-------------------------------|--------------|
| LN <sub>2</sub> Capacity (L)  | 3            |
| Static Hold Time (Days)       | 27           |
| Evap. Rate (Liters/Day)       | 0.11         |
| DIMENSIONS                    |              |
| Necktube (in/mm)              | 2.0 / 50.8   |
| Overall Dia. (in/mm)          | 7.6 / 193    |
| Height (in/mm)                | 17.2 / 436.8 |
| WEIGHT                        |              |
| Empty (lb/kg)                 | 7.2 / 3.2    |
| Full Weight (lb/kg)           | 12.5 / 5.7   |
| CAPACITY                      |              |
| Number of Canisters           | 6            |
| .5 cc Straws on Canes         | N/A          |
| .5 cc Straws – Bulk (1 Level) | 750          |
| Number of 1.2 & 2.0 ml Vials  | N/A          |

# 8XTL

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS                | 8XTL         |
|-------------------------------|--------------|
| LN <sub>2</sub> Capacity (L)  | 8            |
| Static Hold Time (Days)       | 80           |
| Evap. Rate (Liters/Day)       | 0.1          |
| DIMENSIONS                    |              |
| Necktube (in/mm)              | 2.0 / 50.8   |
| Overall Dia. (in/mm)          | 15.6 / 396.2 |
| Height (in/mm)                | 19 / 482.6   |
| WEIGHT                        |              |
| Empty (lb/kg)                 | 19.6 / 8.9   |
| Full Weight (lb/kg)           | 33.8 / 15.3  |
| CAPACITY                      |              |
| Number of Canisters           | 6            |
| .5 cc Straws on Canes         | N/A          |
| .5 cc Straws – Bulk (1 Level) | 750          |
| Number of 1.2 & 2.0 ml Vials  | N/A          |

# 10XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS                | 10XT         |
|-------------------------------|--------------|
| LN <sub>2</sub> Capacity (L)  | 10           |
| Static Hold Time (Days)       | 100          |
| Evap. Rate (Liters/Day)       | 0.1          |
| DIMENSIONS                    |              |
| Necktube (in/mm)              | 2.0 / 50.8   |
| Overall Dia. (in/mm)          | 11.4 / 289.6 |
| Height (in/mm)                | 23.8 / 604.5 |
| WEIGHT                        |              |
| Empty (lb/kg)                 | 16.5 / 7.5   |
| Full Weight (lb/kg)           | 34.3 / 15.6  |
| CAPACITY                      |              |
| Number of Canisters           | 6            |
| .5 cc Straws on Canes         | 540          |
| .5 cc Straws – Bulk (1 Level) | 750          |
| Number of 1.2 & 2.0 ml Vials  | 150          |



# 20XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS                | 20XT         |
|-------------------------------|--------------|
| LN <sub>2</sub> Capacity (L)  | 20.7         |
| Static Hold Time (Days)       | 230          |
| Evap. Rate (Liters/Day)       | 0.09         |
| DIMENSIONS                    |              |
| Necktube (in/mm)              | 2.0 / 50.8   |
| Overall Dia. (in/mm)          | 15.6 / 396.2 |
| Height (in/mm)                | 25.8 / 655.3 |
| WEIGHT                        |              |
| Empty (lb/kg)                 | 26 / 11.8    |
| Full Weight (lb/kg)           | 62.8 / 28.5  |
| CAPACITY                      |              |
| Number of Canisters           | 6            |
| .5 cc Straws on Canes         | 540          |
| .5 cc Straws – Bulk (1 Level) | 750          |
| Number of 1.2 & 2.0 ml Vials  | 150          |

# 34XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS                | 34XT         |
|-------------------------------|--------------|
| LN <sub>2</sub> Capacity (L)  | 34           |
| Static Hold Time (Days)       | 340          |
| Evap. Rate (Liters/Day)       | 0.1          |
| DIMENSIONS                    |              |
| Necktube (in/mm)              | 2.0 / 50.8   |
| Overall Dia. (in/mm)          | 18.8 / 477.5 |
| Height (in/mm)                | 26.3 / 668   |
| WEIGHT                        |              |
| Empty (lb/kg)                 | 34.75 / 15.8 |
| Full Weight (lb/kg)           | 95.3 / 43.2  |
| CAPACITY                      |              |
| Number of Canisters           | 6            |
| .5 cc Straws on Canes         | 540          |
| .5 cc Straws – Bulk (1 Level) | 750          |
| Number of 1.2 & 2.0 ml Vials  | 150          |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.