



designed for scientists



## ULTRA-TURRAX® Tube Drive P control Workstation

/// Data Sheet

The new IKA UTTD power control disperser gets its name from an even more powerful motor. The disperser assures the best stirring, dispersing, grinding and homogenisation results and has been devised in particular for the safe processing of infectious, toxic or odour-intensive samples.

Disposable tubes

The sealable disposable tubes have been optimised for your application.

They reduce the amount of cleaning required and rule out any risk of cross-contamination.



designed for scientists

#### Programmable test conditions

Time, energy, volume: You specify and systematically program the test conditions and store frequently used parameters in your library.

#### Adjustable change in direction of rotation

When, and how often, the direction of rotation should be changed is very easily set using the clear and multilingual menu.

#### Turbo button

Even more intensive mixing, dispersing or grinding at your fingertips: All of this is achieved with the turbo button. Over short periods, rotational speeds of up to 8,000 rpm can be achieved.



designed for scientists

## Technical Data

Motor rating input [W]	36
Motor rating output [W]	28
Viscosity max. [mPas]	5000
Speed range [rpm]	400 - 8000
Reversible direction of rotation	yes
Speed deviation [%]	3
Speed control	10 RPM Steps
Speed display	OLED
Noise without element [dB(A)]	50
Process type	batch
Timer	yes
Timer display	OLED
Dimensions (W x H x D) [mm]	122 x 54 x 178
Weight [kg]	0.9
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 20
USB interface	yes
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input [W]	36
DC Voltage [V=]	24
Current consumption [mA]	1500



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

