


## 96 Well PCR Plate



### PRODUCT SPECIFICATIONS

<b>Catalog Numbers:</b>	IST-401-096MP IST-401-096TP	Pk of 100 Pk of 10	NON-STERILE NON-STERILE
<b>Product Description:</b>	0.2mL 96-well PCR plate Non-skirted, Alpha-numeric labeling, natural-color		
<b>Product Dimensions:</b>	Length: 117.856mm (4.64") Width: 77.724mm (3.06") Height: 20.828mm (.82")		
<b>Packaging:</b>	10 plates per pack 10 packs per case (100 plates total)		
<b>MicroBiological testing:</b>	Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.		
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.		



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.