

**iST-600
384 Well
PCR Plate
extra volume**



PRODUCT SPECIFICATIONS


Catalog Numbers:	<ul style="list-style-type: none"> • IST-600-384MP • IST-600-384TP 	Pk of 100 Pk of 10	NON-STERILE NON-STERILE
Product Description:	The iST Scientific extra volume 384-well full-skirted plate enables post-PCR applications such as enzymatic digestions and ethanol precipitation to be carried out in the wells due to the increased well volume.		
Colour:	Clear		
Product Dimensions:	<ul style="list-style-type: none"> • Length: 127.5mm (5.02") • Width: 85.4mm (3.36") • Height: 10.3mm (0.41") 		
Key Features:	<ul style="list-style-type: none"> • Tall chimney well design • Increased well volume accommodates sequencing and wash steps • Maximum well volume: 40 µL • Usual working well volume: 25 µL • Cut corner: A24 • Polypropylene 		
For Use With (Equipment)	Real-Time PCR Systems or Thermal Cyclers, Thermal Cyclers		
Applications	<ul style="list-style-type: none"> • PCR • qPCR/Real-Time PCR 		
No. of Wells:	384		
Reaction Speed:	Standard		
Skirt Style:	Full		
Volume (Metric) Well:	40 µL		
Volume (Metric) Usual Working:	25 µL		
Type:	PCR Microplate		
Format:	384-well Plate		
Barcode:	Optional		
Packaging:	<ul style="list-style-type: none"> • 10 plates per pack • 10 packs per case (100 plates total) 		
MicroBiological Testing:	Certified testing is performed lot by lot. These tests include Nuclease (Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures		

For Research Use Only. Not for use in diagnostic procedures.

iST-600 384 Well PCR Plate extra volume

Compatibility:

- Agilent Gradient Cycler
- Amersham MegaBace™ 4000
- Biometra T1 Thermocycler
- Bio-Rad C1000
- Bio-Rad CFX384
- Bio-Rad CFX96
- Bio-Rad S1000
- Gene Technologies GS1
- Gene Technologies GS4
- Gene Technologies GSX
- MWG Primus 384
- Techne Flexigene
- Techne Genius
- Techne Prime
- Techne Prime Elite
- Techne PrimeG
- Techne TC-4000
- Techne TC-412
- Techne TC-5000
- Techne TC-512
- Techne TC-PLUS
- Techne Touchgene
- Thermo Scientific Arktik
- Thermo Scientific Hybad MultiBlock System
- Thermo Scientific Hybad PCR Express
- Thermo Scientific Hybad Px2
- Thermo Scientific Hybad PxE





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.